

Datasheet: HCA460

Description:	HUMAN ANTI BUROSUMAB
Specificity:	BUROSUMAB
Other names:	Crysvita
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD63870ia
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG1 antibody (lambda light chain) selected from the HuCAL phage display library and expressed in a human cell line. This antibody is supplied as a liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Source	HKB-11
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml
Immunogen	Burosumab

Specificity **Human Anti-Burosumab Antibody, clone AbD63870ia** is a paratope specific, inhibitory, anti-idiotypic antibody (Type 1) that specifically recognizes the free human monoclonal antibody drug burosumab. It does not recognize the drug target hFGF-23-His, nor burosumab in complex with hFGF-23-His.

HCA460 is a fully human recombinant monoclonal antibody with IgG1 isotype and is suitable as a surrogate positive control or calibrator in an anti-drug antibody (ADA) assay. This antibody has a monovalent intrinsic affinity of 0.6nM and is one of the three antibodies offered for ADA assays. Antibody [HCA459](#) has an affinity of 0.1nM and [HCA461](#) has an affinity of 6.2nM.

Burosumab (trade name Crysvida) is a humanized monoclonal antibody of the IgG1/kappa isotype that binds to fibroblast growth factor-23 (FGF23), blocking the interaction with the FGFR1. Burosumab therefore inhibits FGF23 signaling, thereby increasing tubular phosphate reabsorption and decreasing renal phosphate excretion, as well as increasing serum levels of 1,25(OH)2D and increasing GI absorption of phosphate. As a result, bone mineralization is improved and the risk of bone disease is decreased. Burosumab is indicated for the treatment of X-linked hypophosphatemia and tumor-induced osteomalacia.

[View a summary of all anti-burosumab antibodies](#)

Affinity The monovalent intrinsic affinity of Human anti-Burosumab Antibody, clone AbD63865ia was measured as $K_D = 0.6$ nM by real time, label-free molecular interaction analysis on immobilized burosumab.

ELISA Clone AbD63870ia can be used in indirect ELISA to detect burosumab. It can also be used as a surrogate positive control or a calibrator in an anti-drug antibody (ADA) assay, together with [HCA459](#) and [HCA461](#).

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements This product and/or its use is covered by claims of U.S. patents, and/or pending U.S. and non-U.S. patent applications owned by or under license to Bio-Rad Laboratories, Inc. See bio-rad.com/en-us/trademarks for details.
Crysvida is a trademark of Kyowa Hakko Kirin and Ultragenyx Pharmaceutical.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/HCA460>

Licensed Use For *in vitro*. research purposes and for commercial applications for the provision of in vitro

testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#).

Related Products

Recommended Useful Reagents

[HISPEC ASSAY DILUENT \(BUF049A\)](#)

[HUMAN ANTI BUROSUMAB \(HCA459\)](#)

[HUMAN ANTI BUROSUMAB \(HCA461\)](#)

[HUMAN ANTI BUROSUMAB \(TZA0109\)](#)

[HUMAN ANTI BUROSUMAB \(TZA0110\)](#)

[HUMAN ANTI BUROSUMAB:HRP \(TZA0110P\)](#)

[HUMAN ANTI BUROSUMAB \(DRUG/TARGET COMPLEX\) \(TZA0111\)](#)

[HUMAN ANTI BUROSUMAB \(DRUG/TARGET COMPLEX\):HRP \(TZA0111P\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M450879:260616'

Printed on 16 Jun 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)