

## Datasheet: MCA1133F

**BATCH NUMBER 151321**

<b>Description:</b>	MOUSE ANTI HUMAN CD49a:FITC
<b>Specificity:</b>	CD49a
<b>Other names:</b>	INTEGRIN ALPHA 1 CHAIN, VLA-1
<b>Format:</b>	FITC
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	TS2/7
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/5

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	FITC	490	525
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative</b>	0.09% Sodium Azide		
<b>Stabilisers</b>	1% Bovine Serum Albumin		
<b>Approx. Protein</b>	IgG concentration 0.1 mg/ml		

## Concentrations

---

**Immunogen** Human CTL line directed against HLA-DR.

---

**External Database Links**

**UniProt:**

[P56199](#)    [Related reagents](#)

**Entrez Gene:**

[3672](#) ITGA1    [Related reagents](#)

---

**RRID** AB\_566792

---

**Fusion Partners** Spleen cells from immunized Balb.c mice were fused with cells of the NS1 mouse myeloma.

---

**Specificity**

**Mouse anti Human CD49a antibody, clone TS2/7** recognizes the human alpha 1 integrin sub-unit, which forms the VLA-1 heterodimer in association with the beta 1 integrin. VLA-1 is a receptor for collagen and laminin, and is expressed by chronically activated T cells, melanoma cells and smooth muscle cells.

Mouse anti Human CD49a antibody, clone TS2/7 has been reported as being suitable for use in Western Blotting.

Mouse anti Human CD49a antibody, clone TS2/7 is routinely tested in flow cytometry using HeLa cells.

---

**Flow Cytometry** Use 10ul of the suggested working dilution to label 10<sup>6</sup> cells in 100ul.

---

**References**

1. Hemler, M.E. *et al.* (1984) Glycoproteins of 210,000 and 130,000 m.w. on activated T cells: cell distribution and antigenic relation to components on resting cells and T cell lines. [J Immunol. 132 \(6\): 3011-8.](#)
2. Hemler, M.E. *et al.* (1985) Biochemical characterization of VLA-1 and VLA-2. Cell surface heterodimers on activated T cells. [J Biol Chem. 260 \(28\): 15246-52.](#)
3. Hafler, D.A. *et al.* (1985) Investigation of in vivo activated T cells in multiple sclerosis and inflammatory central nervous system diseases. [Clin Immunol Immunopathol. 37 \(2\): 163-71.](#)
4. Klemmt, P.A. *et al.* (2007) Endometrial cells from women with endometriosis have increased adhesion and proliferative capacity in response to extracellular matrix components: towards a mechanistic model for endometriosis progression. [Hum Reprod. 22: 3139-47.](#)
5. Berry, P.A. *et al.* (2011) The calcium sensor STIM1 is regulated by androgens in prostate stromal cells. [Prostate. 71: 1646-55.](#)
6. Elshaw, S.R. *et al.* (2001) A comparison of ocular melanocyte and uveal melanoma cell invasion and the implication of alpha1beta1, alpha4beta1 and alpha6beta1 integrins. [Br J Ophthalmol. 85: 732-8.](#)
7. Cervar-Zivkovic, M. *et al.* (2011) Endothelin-1 Stimulates Proliferation of First-Trimester Trophoblasts via the A- and B-Type Receptor and Invasion via the B-Type Receptor. [J Clin](#)

[Endocrinol Metab. 96: 3408-15.](#)

8. Kaneko, Y. *et al.* (2012) Integrin  $\alpha 1/\beta 1$  and  $\alpha 2/\beta 1$  as a receptor for IgA1 in human glomerular mesangial cells in IgA nephropathy. [Int Immunol. 24: 219-32.](#)

---

**Storage**

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**

12 months from date of despatch

---

**Health And Safety Information**

Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA1133F>

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA928F\)](#)

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M364825:200529'

**Printed on 12 Mar 2026**