

Datasheet: MCA1283

BATCH NUMBER 173150

Description:	MOUSE ANTI HUMAN CD88
Specificity:	CD88
Other names:	C5aR
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	S5/1
Isotype:	IgG2a
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen		▪		
Immunohistology - Paraffin	▪			
ELISA			▪	
Immunoprecipitation		▪		
Western Blotting	▪			
Functional Assays (1)	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) **Bio-Rad recommend the use of [MCA1283EL](#) for functional studies.**

Target Species

Human

Species Cross Reactivity

Reacts with: Rabbit, Bovine, Ferret, Mink

Based on sequence similarity, is expected to react with: Mustelid

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% sodium azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	C5aR - peptide: Met1 - Asn31.
External Database Links	<p>UniProt: P21730 Related reagents</p> <p>Entrez Gene: 728 C5AR1 Related reagents</p>
Synonyms	C5AR, C5R1
RRID	AB_321883
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the X63-Ag8 myeloma cell line.
Specificity	<p>Mouse anti Human CD88 antibody, clone S5/1 recognizes the C5a receptor (C5aR) CD88, which is predominantly expressed on cells of the myeloid lineage. Clone S5/1 was raised against a synthetic peptide comprising the N-terminal extracellular domain of the C5aR (met1-Asn31) and has recently been shown to recognise the heptameric peptide (D15DKDTLD21).</p> <p>Clone S5/1 has been shown to inhibit the binding of C5a to its receptor.</p>
Flow Cytometry	Use 10µl of the suggested working dilution to label 5 x 10 ⁵ cells in 100µl
References	<ol style="list-style-type: none"> 1. Oppermann, M. & Götze, O. (1994) Plasma clearance of the human C5a anaphylatoxin by binding to leucocyte C5a receptors. Immunology. 82 (4): 516-21. 2. Oppermann, M. <i>et al.</i> (1995) Antibodies from the myeloid panel that react with the C5a receptor and antagonize C5a biological activity. In: Schlossman, S.F. (ed.) Leucocyte Typing V. O.U.P. pp 955-956. 3. Werfel, T. <i>et al.</i> (1996) CD88 antibodies specifically bind to C5aR on dermal CD117+ and CD14+ cells and react with a desmosomal antigen in human skin. J Immunol. 157 (4): 1729-35.

4. Thivierge, M. *et al.* (1999) Modulation of formyl peptide receptor expression by IL-10 in human monocytes and neutrophils. [J Immunol. 162: 3590-5.](#)
5. Eglite, S. *et al.* (2000) Requirements for C5a receptor-mediated IL-4 and IL-13 production and leukotriene C4 generation in human basophils. [J Immunol. 165: 2183-9.](#)
6. Kraft, K. *et al.* (2001) Characterization of sequence determinants within the carboxyl-terminal domain of chemokine receptor CCR5 that regulate signaling and receptor internalization. [J Biol Chem. 276: 34408-18.](#)
7. Sumichika, H. *et al.* (2002) Identification of a potent and orally active non-peptide C5a receptor antagonist. [J Biol Chem. 277: 49403-7.](#)
8. Huang, L. *et al.* (2005) Discovery of human antibodies against the C5aR target using phage display technology. [J Mol Recognit. 18: 327-33.](#)
9. Hüttenrauch, F. *et al.* (2005) G protein-coupled receptor kinases promote phosphorylation and beta-arrestin-mediated internalization of CCR5 homo- and hetero-oligomers. [J Biol Chem. 280: 37503-15.](#)
10. Aasted, B. and Viuff, B. (2007) Reactivity of monoclonal antibodies to human CD antigens with cells from mink. [Vet Immunol Immunopathol. 119: 27-37.](#)
11. Sopp, P. *et al.* (2007) Cross-reactivity of mAbs to human CD antigens with cells from cattle. [Vet Immunol Immunopathol. 119: 106-14.](#)
12. Fukuoka, Y. *et al.* (2008) Generation of anaphylatoxins by human beta-tryptase from C3, C4, and C5. [J Immunol. 180: 6307-16.](#)
13. Schreiber, A. *et al.* (2009) C5a receptor mediates neutrophil activation and ANCA-induced glomerulonephritis. [J Am Soc Nephrol. 20: 289-98.](#)
14. Conroy, A. *et al.* (2009) C5a enhances dysregulated inflammatory and angiogenic responses to malaria in vitro: potential implications for placental malaria. [PLoS One. 4: e4953.](#)
15. Martel, C.J. & Aasted, B. (2009) Characterization of antibodies against ferret immunoglobulins, cytokines and CD markers. [Vet Immunol Immunopathol. 132:109-15.](#)
16. Camous, L. *et al.* (2011) Complement alternative pathway acts as a positive feedback amplification of neutrophil activation. [Blood. 117: 1340-9.](#)
17. Corrales, L. *et al.* (2012) Anaphylatoxin C5a Creates a Favorable Microenvironment for Lung Cancer Progression. [J Immunol. 189: 4674-83.](#)
18. Tseng CW *et al.* (2015) Increased Susceptibility of Humanized NSG Mice to Panton-Valentine Leukocidin and *Staphylococcus aureus* Skin Infection. [PLoS Pathog. 11 \(11\): e1005292.](#)
19. Bettoni, S. *et al.* (2017) Interaction between Multimeric von Willebrand Factor and Complement: A Fresh Look to the Pathophysiology of Microvascular Thrombosis. [J Immunol. 199 \(3\): 1021-40.](#)
20. Tromp, A.T. *et al.* (2020) Host-Receptor Post-Translational Modifications Refine Staphylococcal Leukocidin Cytotoxicity. [Toxins \(Basel\). 12 \(2\): 106.](#)
21. He, J. *et al.* (2023) TPST2-mediated receptor tyrosine sulfation enhances leukocidin cytotoxicity and *S. aureus* infection. [Front Immunol. 14: 1242330.](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in

frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA1283
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M409058:221017'

Printed on 19 Mar 2026