

Datasheet: MCA1289

BATCH NUMBER 174121

Description:	RAT ANTI MOUSE IgG1
Specificity:	IgG1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	LO-MG1-13
Isotype:	IgG1
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			5ug/ml as coating antibody
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	Does not react with:Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.1% sodium azide (NaN ₃)

Stabilisers

Approx. Protein Concentrations	IgG concentration 1 mg/ml
---------------------------------------	---------------------------

Immunogen	Purified mouse IgG1 from BALB/c mice.
------------------	---------------------------------------

External Database Links

UniProt:

[P01869](#) [Related reagents](#)

[P01868](#) [Related reagents](#)

Entrez Gene:

[16017](#) Ighg1 [Related reagents](#)

[16017](#) Ighg1 [Related reagents](#)

Synonyms	Igh-4
-----------------	-------

RRID	AB_321953
-------------	-----------

Fusion Partners	Spleen cells from immunized LOU/c rats were fused with cells of the rat IR983F myeloma cell line.
------------------------	---

Specificity	Rat anti Mouse IgG1 antibody, clone LO-MG1-13 recognizes the gamma 1 heavy chain of mouse immunoglobulin. No cross reactivity with other immunoglobulin classes or subclasses.
--------------------	---

ELISA	This antibody may be used as a coating antibody in a sandwich ELISA in combination with LO-MG1-2 (MCA336P) as the detection reagent.
--------------	--

References

1. Soares, M. *et al.* (1992) Isotype characterization of rat preformed natural antibodies against guinea pig cells. [Transplant Proc. 24 \(2\): 451-2.](#)
 2. Denis, O. *et al.* (1993) Resting B cells can act as antigen presenting cells in vivo and induce antibody responses. [Int Immunol. 5 \(1\): 71-8.](#)
 3. Savignac, M. *et al.* (2010) Increased B cell proliferation and reduced Ig production in DREAM transgenic mice. [J Immunol. 185:7527-36.](#)
 4. Banerjee, K. *et al.* (2012) Occluding the mannose moieties on human immunodeficiency virus type 1 gp120 with griffithsin improves the antibody responses to both proteins in mice. [AIDS Res Hum Retroviruses. 28 \(2\): 206-14.](#)
 5. Laffont, S. *et al.* (2017) Androgen signaling negatively controls group 2 innate lymphoid cells. [J Exp Med. 214 \(6\): 1581-92.](#)
-

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
----------------	---

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10303 available at: https://www.bio-rad-antibodies.com/SDS/MCA1289
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Rat IgG (STAR16...)	DyLight®800
Rabbit Anti Rat IgG (STAR17...)	FITC
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...)	DyLight®550 , DyLight®650 , DyLight®800
Goat Anti Rat IgG (STAR69...)	FITC
Goat Anti Rat IgG (STAR73...)	RPE
Goat Anti Rat IgG (STAR72...)	HRP
Goat Anti Rat IgG (STAR131...)	Alk. Phos. , Biotin
Rabbit Anti Rat IgG (STAR21...)	HRP

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M449783:260213'

Printed on 24 Aug 2026