

## Datasheet: MCA1292

<b>Description:</b>	RAT ANTI MOUSE IgG3
<b>Specificity:</b>	IgG3 HEAVY CHAIN
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	LO-MG3-13
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			5ug/ml (coating)
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on mouse IgG3/Sepharose 4B from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Approx. Protein Concentrations</b>	IgG concentration 1 mg/ml
<b>Immunogen</b>	Polyclonal mouse IgG
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P03987</a> <a href="#">Related reagents</a>

**Entrez Gene:**[380795](#) [AI324046](#) [Related reagents](#)

---

<b>RRID</b>	AB_321848
-------------	-----------

---

<b>Fusion Partners</b>	Spleen cells from immunized Lou/c rats were fused with cells of the rat IR983F myeloma cell line.
------------------------	---

---

<b>Specificity</b>	<b>Rat anti Mouse IgG3 antibody, clone LO-MG3-13</b> recognizes murine IgG3 (gamma heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses.
--------------------	--

---

<b>ELISA</b>	This antibody may be used as a coating antibody in a sandwich ELISA in combination with LO-MG3-7 ( <a href="#">MCA423B</a> ) and Streptavidin:HRP ( <a href="#">STAR5B</a> ) as the detection reagents.
--------------	---

---

<b>References</b>	<ol style="list-style-type: none"><li>1. Denis, O. <i>et al.</i> (1993) Resting B cells can act as antigen presenting cells in vivo and induce antibody responses. <a href="#">Int Immunol. 5 (1): 71-8.</a></li><li>2. Savignac, M. <i>et al.</i> (2010) Increased B cell proliferation and reduced Ig production in DREAM transgenic mice. <a href="#">J Immunol. 185 (12): 7527-36.</a></li><li>3. Banerjee, K. <i>et al.</i> (2009) Enzymatic removal of mannose moieties can increase the immune response to HIV-1 gp120 <i>in vivo</i>. <a href="#">Virology. 389 (1-2): 108-21.</a></li><li>4. Banerjee, K. <i>et al.</i> (2012) Occluding the mannose moieties on human immunodeficiency virus type 1 gp120 with griffithsin improves the antibody responses to both proteins in mice. <a href="#">AIDS Res Hum Retroviruses. 28 (2): 206-14.</a></li></ol>
-------------------	---

---

<b>Storage</b>	Store at +4°C or at -20°C if preferred.  This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: 10040: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a>
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Rat IgG (STAR17...)	<a href="#">FITC</a>
Goat Anti Rat IgG (STAR69...)	<a href="#">FITC</a>
Goat Anti Rat IgG (STAR131...)	<a href="#">Alk. Phos.</a> , <a href="#">Biotin</a>
Goat Anti Rat IgG (STAR73...)	<a href="#">RPE</a>
Rabbit Anti Rat IgG (STAR21...)	<a href="#">HRP</a>
Goat Anti Rat IgG (STAR72...)	<a href="#">HRP</a>
Rabbit Anti Rat IgG (STAR16...)	<a href="#">DyLight®800</a>
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...)	<a href="#">DyLight®800</a>

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

**Printed on 11 Aug 2020**

---