

## Datasheet: MCA1296P

**BATCH NUMBER 151253**

<b>Description:</b>	MOUSE ANTI RAT KAPPA/LAMBDA LIGHT CHAIN:HRP
<b>Specificity:</b>	KAPPA/LAMBDA LIGHT CHAIN
<b>Format:</b>	HRP
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	MARK-1 / MARL-15
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5ug/ml as detecting antibody
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rat
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	50% Glycerol 0.01% Thiomersal
<b>Approx. Protein</b>	IgG concentration 1 mg/ml

## Concentrations

---

**Immunogen** MARK-1: Rat kappa light chains (IR202 and IR968) MARL-15: Rat lambda light chains (RH58, IR31)

---

## External Database

### Links

#### UniProt:

[P01835](#) [Related reagents](#)

[P01836](#) [Related reagents](#)

[P20766](#) [Related reagents](#)

[P20767](#) [Related reagents](#)

#### Entrez Gene:

[500180](#) LOC500180 [Related reagents](#)

[363828](#) RGD1564318 [Related reagents](#)

---

**RRID** AB\_322603

---

**Fusion Partners** Spleen cells from immunized BALB/c mice were fused with cells of the SP2/0 (MARK-1) or PAI-O (MARL-15) mouse myeloma cell line.

---

**Specificity** **Mouse anti Rat Kappa/Lambda Light Chain: antibody cocktail** is a mixture of clones MARK-1 and MARL-15 recognizing rat kappa and lambda light chains respectively. The cocktail does not recognize immunoglobulin heavy chains.

---

**ELISA** HRP-conjugated MARK-1/MARL-15 clone cocktail may be used as a detection reagent for rat immunoglobulins in ELISA assays and can be used in combination with coating antibodies:

rat IgG1	<a href="#">MCA194</a>
rat IgG2a	<a href="#">MCA278</a>
rat IgG2b	<a href="#">MCA1294</a>
rat IgG2c	<a href="#">MCA1295</a>
rat IgM	<a href="#">MCA189</a>
rat IgE	<a href="#">MCA193</a>
rat IgA	<a href="#">MCA191</a>

---

## References

1. Daley, S.R. *et al.* (2007) A key role for TGF-beta signaling to T cells in the long-term acceptance of allografts. [J Immunol. 179 \(6\): 3648-54.](#)
2. Daley, S.R. *et al.* (2008) Fc-disabled anti-mouse CD40L antibodies retain efficacy in promoting transplantation tolerance. [Am J Transplant. 8 \(11\): 2265-71.](#)
3. Sun, Q. *et al.* (2012) Upregulated protein arginine methyltransferase 1 by IL-4 increases eotaxin-1 expression in airway epithelial cells and participates in antigen-induced pulmonary inflammation in rats. [J Immunol. 188 \(7\): 3506-12.](#)
4. Cousins, L. *et al.* (2006) Eosinophilic bowel disease controlled by the BB rat-derived lymphopenia/Gimap5 gene. [Gastroenterology. 131:1475-85.](#)
5. Sandin, L.C. *et al.* (2014) Locally delivered CD40 agonist antibody accumulates in

secondary lymphoid organs and eradicates experimental disseminated bladder cancer.  
[Cancer Immunol Res. 2 \(1\): 80-90.](#)

---

**Storage** Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10097 available at:  
<https://www.bio-rad-antibodies.com/SDS/MCA1296P>  
10097

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M365103:200529'

Printed on 25 Mar 2023

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)