

Datasheet: MCA1366GA

BATCH NUMBER 158468

Description:	HAMSTER ANTI MOUSE TCR GAMMA/DELTA
Specificity:	TCR GAMMA/DELTA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	GL3
Isotype:	IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/5 - 1/50
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 0.5mg/ml

Immunogen	C57BL/6J intra epithelial lymphocytes.
RRID	AB_1102884
Fusion Partners	Spleen cells from immunized Armenian hamster were fused with cells of the murine myeloma SP2/0 myeloma cell line.
Specificity	<p>Hamster anti Mouse TCR Gamma/Delta antibody, clone GL3 reacts with mouse TCR gamma/delta expressing lymphocytes.</p> <p>Hamster anti Mouse TCR Gamma/Delta antibody, clone GL3 has been shown to have depleting activity <i>in vivo</i> (Van der Heyde et al. 1995 and Skeen et al. 1993).</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> 1. Skarstein, K. <i>et al.</i> (1994) Oligoclonality of T cells in salivary glands of autoimmune MRL/lpr mice. Immunology. 81 (4): 497-501. 2. van der Heyde, H.C. <i>et al.</i> (1995) Gamma delta T cells function in cell-mediated immunity to acute blood-stage Plasmodium chabaudi adami malaria. J Immunol. 154 (8): 3985-90. 3. Skeen, M.J. & Ziegler, H.K. (1993) Induction of murine peritoneal gamma/delta T cells and their role in resistance to bacterial infection. J Exp Med. 178 (3): 971-84. 4. Karasova, D. <i>et al.</i> (2010) Influence of 5 major Salmonella pathogenicity islands on NK cell depletion in mice infected with <i>Salmonella enterica</i> serovar Enteritidis. BMC Microbiol. 10: 75. 5. Laky, K. <i>et al.</i> (2000) Enterocyte expression of interleukin 7 induces development of gammadelta T cells and Peyer's patches. J Exp Med. 191: 1569-80.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>
Guarantee	Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA1366GA 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Hamster IgG (STAR104...) [DyLight®550](#), [DyLight®650](#), [DyLight®800](#),

FITC

Goat Anti Hamster IgG (STAR79...) Biotin, FITC, HRP

Recommended Negative Controls

HAMSTER (ARMENIAN) IgG NEGATIVE CONTROL (MCA2356)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M382919:210513'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)