

## Datasheet: MCA1400G

**BATCH NUMBER 0512**

<b>Description:</b>	MOUSE ANTI HUMAN FGF BASIC
<b>Specificity:</b>	FGF BASIC
<b>Other names:</b>	FGF2
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	MC-GF1
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.5 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide

<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml
<b>Immunogen</b>	Recombinant human bFGF.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P09038</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">2247</a>    FGF2    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	FGFB
<b>RRID</b>	AB_2278442
<b>Fusion Partners</b>	Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2 myeloma cell line.
<b>Specificity</b>	<p><b>Mouse anti Human FGF Basic antibody, clone MC-GF1</b> recognizes basic fibroblast growth factor (bFGF). Mouse anti Human FGF Basic antibody, clone MC-GF1 binds strongly to bFGF attached to proteoglycans.</p> <p>Mouse anti Human FGF Basic antibody, clone MC-GF1 is not suitable for detecting bFGF bound to cell surface receptors.</p>
<b>References</b>	<ol style="list-style-type: none"> <li>Ribatti, D. <i>et al.</i> (2003) Neovascularisation, expression of fibroblast growth factor-2, and mast cells with tryptase activity increase simultaneously with pathological progression in human malignant melanoma. <a href="#">Eur J Cancer. 39:666-74.</a></li> <li>Whitworth, M.K. <i>et al.</i> (2005) Regulation of fibroblast growth factor-2 activity by human ovarian cancer tumor endothelium. <a href="#">Clin Cancer 11: 4282-8.</a></li> <li>Zong, F. <i>et al.</i> (2009) Syndecan-1 and FGF-2, but not FGF receptor-1, share a common transport route and co-localize with heparanase in the nuclei of mesenchymal tumor cells. <a href="#">PLoS One. 4:e7346.</a></li> <li>Melrose, J. <i>et al.</i> (2012) Chondroitin sulphate and heparan sulphate sulphation motifs and their proteoglycans are involved in articular cartilage formation during human foetal knee joint development. <a href="#">Histochem Cell Biol. 138 (3): 461-75.</a></li> </ol>
<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch

---

**Health And Safety Information**      Material Safety Datasheet documentation #10040 available at:  
<https://www.bio-rad-antibodies.com/SDS/MCA1400G10040>

---

**Regulatory**                      For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)                      [HRP](#)  
Rabbit Anti Mouse IgG (STAR12...)                      [RPE](#)  
Goat Anti Mouse IgG (STAR70...)                      [FITC](#)  
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)  
Goat Anti Mouse IgG (STAR76...)                      [RPE](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR13...)                      [HRP](#)  
Goat Anti Mouse IgG (Fc) (STAR120...)                      [FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...)                      [FITC](#)

**North & South America**      Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M365268:200529'

**Printed on 25 Mar 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)