

## Datasheet: MCA152P

**BATCH NUMBER 158586**

<b>Description:</b>	RAT ANTI MOUSE KAPPA LIGHT CHAIN:HRP
<b>Specificity:</b>	KAPPA LIGHT CHAIN
<b>Format:</b>	HRP
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	OX-20
<b>Isotype:</b>	IgG2b
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/100 - 1/10,000
Immunoprecipitation			▪	
Western Blotting	▪			1/1000 - 1/10,000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal HRP Stabiliser ( <a href="#">BUF052A</a> )
<b>Approx. Protein</b>	IgG concentration 1.0 mg/ml

## Concentrations

---

**Immunogen** Purified mouse IgG.

---

## External Database Links

**UniProt:**

[P01837](#)    [Related reagents](#)

**Entrez Gene:**

[16071](#) Igk-C    [Related reagents](#)

---

**Fusion Partners** Spleen cells from an immunized PVG.RT1c rat were fused with cells of the rat Y3 Ag.1.2.3 myeloma cell line.

---

**Specificity** **Rat anti Mouse kappa light chain antibody, clone OX-20** recognizes the kappa light chain of murine immunoglobulins. OX-20 does not recognize mouse lambda light chains or bind to the immunoglobulin heavy chain.

---

## References

1. Mason, D.W. *et al.* (1981) The rat mixed lymphocyte reaction: roles of a dendritic cell in intestinal lymph and T-cell subsets defined by monoclonal antibodies. [Immunology. 44 \(1\): 75-87.](#)
  2. White, R.A. *et al.* (1984) Differentiation markers in the Drosophila ovary. [J Embryol Exp Morphol. 84: 275-86.](#)
  3. Meyer, K.B. *et al.* (1990) The importance of the 3'-enhancer region in immunoglobulin kappa gene expression. [Nucleic Acids Res. 18 \(19\): 5609-15.](#)
  4. Thomas, M.J. *et al.* (1999) CD8 T-lymphocyte-mediated regulation of ovalbumin-specific murine IgE responses. [Int Arch Allergy Immunol. 118 \(2-4\): 289-91.](#)
  5. Wasserman, G.M. *et al.* (1993) An approach to development of specific T-lymphocyte lines by use of preprocessed antigens in *Plasmodium vinckei vinckei* murine malaria. [Infect Immun. 61 \(5\): 1958-63.](#)
  6. Thomas, M.J. *et al.* (2002) CD8 T cells inhibit IgE via dendritic cell IL-12 induction that promotes Th1 T cell counter-regulation. [J Immunol. 168 \(1\): 216-23.](#)
  7. Salagianni, M. *et al.* (2007) An essential role for IL-18 in CD8 T cell-mediated suppression of IgE responses. [J Immunol. 178 \(8\): 4771-8.](#)
- 

**Further Reading** 1. Galfrè, G. *et al.* (1979) Rat x rat hybrid myelomas and a monoclonal anti-Fd portion of mouse IgG. [Nature. 277 \(5692\): 131-3.](#)

---

## Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

---

## Guarantee

12 months from date of despatch

---

## Health And Safety Information

Material Safety Datasheet documentation #10131 available at: <https://www.bio-rad-antibodies.com/SDS/MCA152P>  
10131

---

## Related Products

### Recommended Useful Reagents

[TidyBlot WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M365386:200529'

**Printed on 18 Jan 2024**

---