

Datasheet: MCA1553F

BATCH NUMBER 174436

Description:	RAT ANTI HUMAN INTERLEUKIN-2:FITC
Specificity:	IL-2
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	MQ1-17H12
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Membrane permeabilization is required for this application. The use of Leucoperm (Product Code [BUF09](#)) is recommended for this purpose.

Target Species	Human								
Species Cross Reactivity	Reacts with: Cynomolgus monkey N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.								
Product Form	Purified IgG - liquid								
Max Ex/Em	<table><tr><th>Fluorophore</th><th>Excitation Max (nm)</th><th>Emission Max (nm)</th></tr><tr><td>FITC</td><td>490</td><td>525</td></tr></table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	FITC	490	525		
Fluorophore	Excitation Max (nm)	Emission Max (nm)							
FITC	490	525							
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant								

Buffer Solution	Phosphate buffered saline
Preservative	0.09% sodium azide (NaN ₃)
Stabilisers	1% bovine serum albumin
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
External Database Links	UniProt: P60568 Related reagents Entrez Gene: 3558 IL2 Related reagents
RRID	AB_2125584
Specificity	Rat anti Human Interleukin-2 antibody, clone MQ1-17H12 recognizes human interleukin-2 (IL-2), also known as T-cell growth factor or Aldesleukin. IL-2 is a 153 amino acid secreted cytokine produced by T-cells following antigenic or mitogenic stimulation. It plays a crucial role in T-cell proliferation and immune regulation. Rat anti Human Interleukin-2 antibody, clone MQ1-17H12 has been successfully used for the detection and measurement of native and recombinant human IL-2 using flow cytometry and ELISA where it acts as the capture reagent in a very sensitive assay for IL-2 determination in human samples. IL-2 is used extensively for the treatment of patients with metastatic melanoma (Atkins et al. 1999) and renal cancer with reproducible beneficial effects (Rosenberg 2014). Monitoring of IL-2 levels in such patients is required for optimization of treatment regimens (Lotze et al. 1986). Rat anti Human Interleukin-2 antibody, clone MQ1-17H12 is also useful for evaluating levels of IL-2 in patients with autoimmune conditions including lupus erythematosus and rheumatoid arthritis where there may be an impaired ability to produce IL-2 (Talal et al. 1983).
Flow Cytometry	Use 10µl of the suggested working dilution to label 1x10 ⁶ cells in 100µl
References	1. Garbe, Y. <i>et al.</i> (2011) Semiallogenic fusions of MSI(+) tumor cells and activated B cells induce MSI-specific T cell responses. BMC Cancer. 11: 410. 2. Moncrieffe, H. <i>et al.</i> (2010) High expression of the ectonucleotidase CD39 on T cells from the inflamed site identifies two distinct populations, one regulatory and one memory T cell population. J Immunol. 185 (1): 134-43.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10041 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA1553F>

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL:FITC \(MCA6005F\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M410032:221024'

Printed on 24 Aug 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)