

Datasheet: MCA1585

BATCH NUMBER 174198

Description:	MOUSE ANTI HUMAN INTEGRIN ALPHA 9 BETA 1
Specificity:	INTEGRIN ALPHA 9 BETA 1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	Y9A2
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen		▪		
Immunohistology - Paraffin		▪		
Immunohistology - Resin		▪		
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting		▪		

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Reacts with: Pig, Guinea Pig N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture

supernatant

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% sodium azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	alpha9 integrin transfected lambda cells.
External Database Links	UniProt: Q13797 Related reagents Entrez Gene: 3680 ITGA9 Related reagents
RRID	AB_322166
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Specificity	Mouse anti Human Integrin alpha 9 beta 1 antibody, clone Y9A2 recognizes the human alpha9 integrin, a ~140 kDa glycoprotein expressed as a heterodimer in association with the beta1 integrin. The alpha9 integrin is expressed by epithelial cells, including airway epithelial cells. Mouse anti Human Integrin alpha 9 beta 1 antibody, clone Y9A2 inhibits adhesion of alpha9 expressing cells to N-terminal osteopontin. Mouse anti Human Integrin alpha 9 beta 1 antibody, clone Y9A2 has been shown to cause 100% inhibition of α9β1-mediated adhesion to the third fibronectin type III repeat in tenascin at <10ug/ml (Yokosaki <i>et al.</i> 1998) and demonstrates Itga9 blocking of progenitor cell proliferation (Schreiber <i>et al.</i> 2009).
Flow Cytometry	Use 10µl of the suggested working dilution to label 10 ⁶ cells in 100µl
References	1. Wang, A. <i>et al.</i> (1996) Differential regulation of airway epithelial integrins by growth factors. Am J Respir Cell Mol Biol. 15 (5): 664-72. 2. Smith, L.L. <i>et al.</i> (1996) Osteopontin N-terminal domain contains a cryptic adhesive sequence recognized by alpha9beta1 integrin. J Biol Chem. 271 (45): 28485-91. 3. Yokosaki, Y. <i>et al.</i> (1996) Differential effects of the integrins alpha9beta1, alphavbeta3, and alphavbeta6 on cell proliferative responses to tenascin. Roles of the beta subunit extracellular and cytoplasmic domains. J Biol Chem. 271 (39): 24144-50. 4. Sime, W. <i>et al.</i> (2009) Human mast cells adhere to and migrate on epithelial and vascular basement membrane laminins LM-332 and LM-511 via alpha3beta1 integrin. J

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA1585
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M437898:250320'

Printed on 24 Aug 2026