

Datasheet: MCA1588P

BATCH NUMBER 173220

Description:	RAT ANTI MOUSE IgG2a HEAVY CHAIN:HRP
Specificity:	IgG2a HEAVY CHAIN
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	LO-MG2a-9
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			250ng/ml
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG conjugated to horseradish peroxidase (HRP) - liquid
Preparation	Purified IgG prepared by affinity chromatography from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.0095% MIT 50% Glycerol
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Purified IgG2a from Balb/c mice. PTR

External Database**Links****UniProt:**

[P01865](#) [Related reagents](#)
[P01864](#) [Related reagents](#)
[P01863](#) [Related reagents](#)

Entrez Gene:

[380793](#) Igh-1a [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)

RRID

AB_322872

Fusion Partners

Spleen cells from an immunized LOU/c rat were fused with cells of the IR983F rat myeloma cell line.

Specificity

Rat anti Mouse IgG2a Heavy Chain antibody, clone LO-MG2a-9 recognizes the mouse gamma 2a heavy chain, binding to all allotypes and does not cross-react with other murine immunoglobulin classes or subclasses.

The avidity of Rat anti Mouse IgG2a Heavy Chain antibody, clone LO-MG2a-9 for murine IgG2a = $2.5 \times 10^9 \text{M}^{-1}$.

References

1. Echeverria, P.C. *et al.* (2006) Potent antigen-specific immunity to *Toxoplasma gondii* in adjuvant-free vaccination system using Rop2-*Leishmania infantum* Hsp83 fusion protein. [Vaccine. 24: 4102-10.](#)
2. Hajek, A.R. *et al.* (2008) 12/15-Lipoxygenase deficiency protects mice from allergic airways inflammation and increases secretory IgA levels. [J Allergy Clin Immunol. 122: 633-9.e3.](#)
3. Savignac, M. *et al.* (2010) Increased B cell proliferation and reduced Ig production in DREAM transgenic mice. [J Immunol. 185:7527-36.](#)
4. Blodörn, K. *et al.* (2014) Vaccine safety and efficacy evaluation of a recombinant bovine respiratory syncytial virus (BRSV) with deletion of the SH gene and subunit vaccines based on recombinant human RSV proteins: N-nanorings, P and M2-1, in calves with maternal antibodies. [PLoS One. 9 \(6\): e100392.](#)
5. Carew, J.S. *et al.* (2015) Targeting Survivin Inhibits Renal Cell Carcinoma Progression and Enhances the Activity of Temsirolimus. [Mol Cancer Ther. 14 \(6\): 1404-13.](#)
6. Islam, S. *et al.* (2020) Resistance to histone deacetylase inhibitors confers hypersensitivity to oncolytic reovirus therapy. [Blood Adv. 4 \(20\): 5297-310.](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #20516 available at: https://www.bio-rad-antibodies.com/SDS/MCA1588P
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M444652:250714'

Printed on 31 Mar 2026