

Datasheet: MCA1658

BATCH NUMBER 166426

Description:	MOUSE ANTI SHEEP INTERLEUKIN-1 BETA
Specificity:	IL-1 BETA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	1D4
Isotype:	IgG1
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	▪			1/10
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Membrane permeabilization is required for this application. The use of Leucoperm (Product Code [BUF09](#)) is recommended for this purpose.

Target Species	Sheep
Species Cross Reactivity	<p>Reacts with: Goat, Bovine, Horse</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid

Preparation	Antibody purified from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	<0.1% sodium azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant ovine IL-1 Beta
External Database Links	<p>UniProt: P21621 Related reagents</p> <p>Entrez Gene: 443539 IL-1B Related reagents</p>
RRID	AB_322126
Specificity	<p>Mouse anti Sheep Interleukin-1 beta antibody, clone 1D4 recognizes ovine interleukin-1 beta, and shows no cross-reactivity with ovine IL-6, IL-8, TNF alpha or MCP-1.</p> <p>Mouse anti Sheep Interleukin-1 beta antibody, clone 1D4 demonstrates partial neutralizing activity of ovine IL-1 beta.</p>
Flow Cytometry	Use 10µl of the suggested working dilution to label 1x10 ⁶ cells in 100µl
ELISA	<p>Mouse anti interleukin-1 beta antibody, clone 1D4 may be used as a capture antibody in a bovine IL-1 beta sandwich ELISA together with Rabbit anti Bovine interleukin-1β antibody (AHP851B) as the detection reagent for evaluation of IL-1β levels in bovine samples together with recombinant Bovine interleukin-1β (PBP008) used as standards.</p> <p>Alternatively, Mouse anti Interleukin-1 beta antibody, clone 1D4 can be used as a capture reagent together with Rabbit anti Ovine interleukin-1β antibody (AHP423) as a detection reagent for the evaluation of IL-1β levels in ovine, bovine or caprine samples, again utilizing recombinant bovine IL-1&beta (PBP008) as an internal standard.</p>
References	<ol style="list-style-type: none"> 1. Martoriati, A. & Gérard, N. (2003) Interleukin-1 (IL-1) system gene expression in granulosa cells: kinetics during terminal preovulatory follicle maturation in the mare. Reprod Biol Endocrinol. 1: 42-51. 2. Bannerman, D.D. <i>et al.</i> (2004) <i>Escherichia coli</i> and <i>Staphylococcus aureus</i> elicit differential innate immune responses following intramammary infection. Clin Diagn Lab Immunol. 11: 463-72. 3. Bannerman, D.D. <i>et al.</i> (2004) Characterization of the bovine innate immune response to intramammary infection with <i>Klebsiella pneumoniae</i>. J Dairy Sci. 87 (8): 2420-32. 4. Leite, F. <i>et al.</i> (2005) Incubation of bovine PMNs with conditioned medium from BHV-1 infected peripheral blood mononuclear cells increases their susceptibility to <i>Mannheimia</i>

- haemolytica* leukotoxin. [Vet Immunol Immunopathol. 103 \(3-4\): 187-93.](#)
5. Matthews, K. *et al.* (2007) The effect of gene gun-delivered pGM-CSF on the immunopathology of the vaccinated skin. [Scand J Immunol. 65 \(3\): 298-307.](#)
 6. Jacobsen, S. *et al.* (2007) The cytokine response of circulating peripheral blood mononuclear cells is changed after intravenous injection of lipopolysaccharide in cattle. [Vet J. 174 \(1\): 170-5.](#)
 7. Cox, R.A. *et al.* (2007) Production of pro-inflammatory polypeptides by airway mucous glands and its potential significance. [Pulm Pharmacol Ther. 20 \(2\): 172-7.](#)
 8. Rainard P *et al.* (2008) *Staphylococcus aureus* lipoteichoic acid triggers inflammation in the lactating bovine mammary gland. [Vet Res. 39 \(5\): 52.](#)
 9. Bougarn, S. *et al.* (2010) Muramyl dipeptide synergizes with *Staphylococcus aureus* lipoteichoic acid to recruit neutrophils in the mammary gland and to stimulate mammary epithelial cells. [Clin Vaccine Immunol. 17 \(11\): 1797-809.](#)
 10. Rinaldi, M. *et al.* (2010) A sentinel function for teat tissues in dairy cows: dominant innate immune response elements define early response to *E. coli* mastitis. [Funct Integr Genomics. 10: 21-38.](#)
 11. Simojoki, H. *et al.* (2011) Innate immune response in experimentally induced bovine intramammary infection with *Staphylococcus simulans* and *S. epidermidis*. [Vet Res. 42: 49.](#)
 12. Redondo, E. *et al.* (2014) Induction of interleukin-8 and interleukin-12 in neonatal ovine lung following experimental inoculation of bovine respiratory syncytial virus. [J Comp Pathol. 150 \(4\): 434-48.](#)
 13. Karisnani K *et al.* (2015) Interleukin-1 Receptor Antagonist Protects against Lipopolysaccharide Induced Diaphragm Weakness in Preterm Lambs. [PLoS One. 10 \(4\): e0124390.](#)
 14. Doull, L. *et al.* (2015) Late production of CXCL8 in ruminant oro-nasal turbinate cells in response to *Chlamydia abortus* infection. [Vet Immunol Immunopathol. 168 \(1-2\): 97-102.](#)
 15. Xu, A. *et al.* (2015) The Ovine Fetal and Placental Inflammatory Response to Umbilical Cord Occlusions With Worsening Acidosis. [Reprod Sci. 22 \(11\): 1409-20.](#)
 16. Sobotta, K. *et al.* (2016) *Coxiella burnetii*. Infects Primary Bovine Macrophages and Limits Their Host Cell Response. [Infect Immun. 84 \(6\): 1722-34.](#)
 17. Cortes, M. *et al.* (2017) RNAseq profiling of primary microglia and astrocyte cultures in near-term ovine fetus: A glial *in vivo-in vitro* multi-hit paradigm in large mammalian brain. [J Neurosci Methods. 276: 23-32.](#)
 18. Canal AM *et al.* (2017) Immunohistochemical detection of pro-inflammatory and anti-inflammatory cytokines in granulomas in cattle with natural *Mycobacterium bovis* infection. [Res Vet Sci. 110: 34-39.](#)
 19. Cao, M. *et al.* (2019) $\alpha 7$ Nicotinic Acetylcholine Receptor Signaling Modulates Ovine Fetal Brain Astrocytes Transcriptome in Response to Endotoxin. [Front Immunol. 10: 1063.](#)
 20. Stassi, A.F. *et al.* (2019) Follicular structures of cows with cystic ovarian disease present altered expression of cytokines. [Zygote. 15: 1-14.](#)
 21. Ciliberti, M.G. *et al.* (2022) Green extraction of bioactive compounds from wine lees and their bio-responses on immune modulation using *in vitro*. sheep model. [J Dairy Sci. Mar 17 \[Epub ahead of print\].](#)
 22. Damiano, S. *et al.* (2022) Red orange and lemon extract preserve from oxidative stress, DNA damage and inflammatory status in lambs [It J Anim Sci. 21 \(1\): 934-42.](#)
 23. Stassi, A.F. *et al.* (2018) Altered expression of IL-1 β , IL-1RI, IL-1RII, IL-1RA and IL-4

could contribute to anovulation and follicular persistence in cattle. [Theriogenology. 110: 61-73.](#)

24. Ciliberti, M.G. *et al.* (2024) Role of hazelnut skin supplementation on plasma antioxidant status and cytokine profile in growing lambs. [Front Vet Sci. 11: 1340141.](#)

Further Reading	1. Rothel, J.S. <i>et al.</i> (1997) Analysis of ovine IL-1 beta production <i>in vivo</i> and <i>in vitro</i> by enzyme immunoassay and immunohistochemistry. Vet Immunol Immunopathol. 57 (3-4): 267-78.
------------------------	--

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
----------------	---

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA1658 10040
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)

Goat Anti Mouse IgG (H/L) (STAR117...) [FITC](#)

Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M421731:230714'

Printed on 04 Mar 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)