

Datasheet: MCA1662

BATCH NUMBER 170081

Description:	MOUSE ANTI HUMAN IgM
Specificity:	IgM
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	M15/8
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			20ug/ml
Immunohistology - Frozen	▪			20ug/ml
Immunohistology - Paraffin			▪	
ELISA	▪			5ug/ml as coating Ab
Immunoprecipitation	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	TRIS buffered saline.
Preservative Stabilisers	<0.1% sodium azide (NaN ₃)
Approx. Protein	IgG concentration 1.0 mg/ml

Concentrations

External Database

Links

UniProt:

[P01871](#)

[Related reagents](#)

Entrez Gene:

[3507](#)

IGHM

[Related reagents](#)

RRID

AB_321861

Fusion Partners

Spleen cells from immunized BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.

Specificity

Mouse anti Human IgM monoclonal antibody, clone M15/8 recognizes human IgM, binding to the μ heavy chain. The antibody does not cross-react with other immunoglobulin classes.

References

1. Cragg, M.S. *et al.* (1999) Analysis of the interaction of monoclonal antibodies with surface IgM on neoplastic B-cells. [Br J Cancer. 79 \(5-6\): 850-7.](#)
2. Irie, E. *et al.* (2010) Severe hypogammaglobulinemia persisting for 6 years after treatment with rituximab combined chemotherapy due to arrest of B lymphocyte differentiation together with alteration of T lymphocyte homeostasis. [Int J Hematol. 91: 501-8.](#)
3. Creasey, A.M. *et al.* (2003) Nonspecific immunoglobulin M binding and chondroitin sulfate A binding are linked phenotypes of *Plasmodium falciparum* isolates implicated in malaria during pregnancy. [Infect Immun. 71: 4767-71.](#)
4. Ghumra, A. *et al.* (2012) Induction of strain-transcending antibodies against Group A PfEMP1 surface antigens from virulent malaria parasites. [PLoS Pathog. 8 \(4\): e1002665.](#)
5. Semblat, J.P. *et al.* (2015) Identification of the minimal binding region of a Plasmodium falciparum IgM binding PfEMP1 domain. [Mol Biochem Parasitol. 201 \(1\): 76-82.](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10057 available at: <https://www.bio-rad-antibodies.com/SDS/MCA1662>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M410983:221031'

Printed on 01 Apr 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)