

Datasheet: MCA1730F

Description:	MOUSE ANTI HUMAN CD44v6:FITC
Specificity:	CD44v6
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	VFF-7
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Buffer Solution	Tris buffered saline		
Preservative	0.02% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.1mg/ml		
Immunogen	Glutathione S Transferase (GST) fusion protein corresponding to the variable domains (v3		

to v10) of human CD44.

External Database**Links****UniProt:**[P16070](#)[Related reagents](#)**Entrez Gene:**[960](#)

CD44

[Related reagents](#)

Synonyms

LHR, MDU2, MDU3, MIC4

RRID

AB_322694

Fusion Partners

Spleen cells from immunized BALB/c micewere fused with cells of the P3X63Ag8.653 myeloma cell line.

Specificity

Mouse anti Human CD44v6 antibody, clone VFF-7 recognizes an epitope encoded by exon v6 on the variant portion of human CD44.

Flow Cytometry

Use 10ul of the suggested working dilution to label 10^6 cells in 100ul.

References

1. Neumayer, R. *et al.* (1999) CD44 expression in benign and malignant colorectal polyps. [Dis Colon Rectum. 42 \(1\): 50-5.](#)
2. Koopman, G. *et al.* (1993) Activated human lymphocytes and aggressive non-Hodgkin's lymphomas express a homologue of the rat metastasis-associated variant of CD44. [J Exp Med. 177 \(4\): 897-904.](#)
3. Yoo, S.H. *et al.* (2005) Expression patterns of markers for type II pneumocytes in pulmonary sclerosing hemangiomas and fetal lung tissues. [Arch Pathol Lab Med. 129 \(7\): 915-9.](#)
4. Kim, J.S. and Werth, V.P. (2011) Identification of specific chondroitin sulfate species in cutaneous autoimmune disease. [J Histochem Cytochem. 59: 780-90.](#)
5. Chaiyarit, P. *et al.* (2008) Alteration of the expression of CD44 isoforms in oral epithelia and saliva from patients with oral lichen planus. [J Clin Immunol. 28: 26-34.](#)
6. Shirure, V.S. *et al.* (2015) CD44 variant isoforms expressed by breast cancer cells are functional E-selectin ligands under flow conditions. [Am J Physiol Cell Physiol. 308 \(1\): C68-78.](#)

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10058 available at:
10058: <https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf>

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA928F\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M365632:200529'

Printed on 12 Feb 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)