

## Datasheet: MCA1734

**BATCH NUMBER 153622**

<b>Description:</b>	MOUSE ANTI HUMAN CD44v3
<b>Specificity:</b>	CD44v3
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	VFF-327
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			10ug/ml
Immunohistology - Frozen	▪			10ug/ml
Immunohistology - Paraffin (1)	▪			10ug/ml
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.**

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline

<b>Preservative Stabilisers</b>	0.09% Sodium Azide
<b>Approx. Protein Concentrations</b>	IgG concentration 1 mg/ml
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P16070</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">960</a>    CD44    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	LHR, MDU2, MDU3, MIC4
<b>RRID</b>	AB_322700
<b>Specificity</b>	<b>Mouse anti Human CD44v3 antibody, clone VFF-327</b> recognizes an epitope encoded by exon V3 on the variant portion of human CD44.
<b>Flow Cytometry</b>	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells or 100ul whole blood
<b>References</b>	<ol style="list-style-type: none"> <li>1. Alam, T.N. <i>et al.</i> (2004) Differential expression of CD44 during human prostate epithelial cell differentiation. <a href="#">J Histochem Cytochem. 52 (8): 1083-90.</a></li> <li>2. Kim, J.S. &amp; Werth, V.P. (2011) Identification of specific chondroitin sulfate species in cutaneous autoimmune disease. <a href="#">J Histochem Cytochem. 59 (8): 780-90.</a></li> <li>3. Bauer, S. <i>et al.</i> (2006) Fibroblast activation protein is expressed by rheumatoid myofibroblast-like synoviocytes. <a href="#">Arthritis Res Ther. 8: R171.</a></li> <li>4. Hanley, W.D. <i>et al.</i> (2006) Variant isoforms of CD44 are P- and L-selectin ligands on colon carcinoma cells. <a href="#">FASEB J. 20: 337-9.</a></li> <li>5. Shirure, V.S. <i>et al.</i> (2015) CD44 variant isoforms expressed by breast cancer cells are functional E-selectin ligands under flow conditions. <a href="#">Am J Physiol Cell Physiol. 308 (1): C68-78.</a></li> <li>6. Spadea, A. <i>et al.</i> (2019) Evaluating the Efficiency of Hyaluronic Acid for Tumor Targeting via CD44. <a href="#">Mol Pharm. 16 (6): 2481-93.</a></li> <li>7. Birzele, F. <i>et al.</i> (2015) Markers for the responsiveness to anti-cd44 antibodies <a href="#">Patent Application Publication: US 2015/0044232 A1</a></li> </ol>
<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety</b>	Material Safety Datasheet documentation #10040 available at:

**Information** <https://www.bio-rad-antibodies.com/SDS/MCA1734>  
10040

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)  
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)  
Goat Anti Mouse IgG (STAR76...) [RPE](#)  
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)  
Goat Anti Mouse IgG (STAR70...) [FITC](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),  
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),  
[FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)  
Goat Anti Mouse IgG (STAR77...) [HRP](#)  
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M365635:200529'

Printed on 18 Jan 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)