

## Datasheet: MCA1744F

**BATCH NUMBER 159377**

<b>Description:</b>	MOUSE ANTI HUMAN CD66e:FITC
<b>Specificity:</b>	CD66e
<b>Other names:</b>	CEA
<b>Format:</b>	FITC
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	C365D3 (NCRC23)
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	FITC	490	525
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative</b>	0.09% Sodium Azide		
<b>Stabilisers</b>	1% Bovine Serum Albumin		
<b>Approx. Protein</b>	IgG concentration 0.1 mg/ml		

## Concentrations

---

### External Database

#### Links

#### UniProt:

[P06731](#)    [Related reagents](#)

#### Entrez Gene:

[1048](#)    CEACAM5    [Related reagents](#)

---

#### Synonyms

CEA

---

#### RRID

AB\_323914

---

#### Fusion Partners

Spleen cells from immunised BALB/c mice were fused with cells of the mouse P3NSI myeloma cell line.

---

#### Specificity

**Mouse anti Human CD66e antibody, clone C365D3 (NCRC23)** recognizes human Carcinoembryonic antigen-related cell adhesion molecule 5, also known as CD66e, carcinoembryonic antigen, Meconium antigen 100, CEA or CEACAM5. CD66e is a 702 amino acid ~77 kDa GPI anchored membrane protein containing 7 [Ig-like domains](#). Mouse anti Human CD66e antibody, clone C365D3 does not cross-react with normal cross-reacting antigen (CD66c), or with biliary glycoprotein 1 (CD66a) as indicated by binding assays ([Price 1988](#), note: in this study Mouse anti Human CD66e antibody, clone C365D3 is designated as clone 6 (from author)).

---

#### Flow Cytometry

Use 10ul of the suggested working dilution to label 10<sup>6</sup> cells in 100ul.

---

#### References

1. Seth, J. *et al.* (1988) Carcinoembryonic antigen. [Lancet. 1 \(8599\): 1399.](#)
2. Soucek, K. *et al.* (2010) Fetal colon cell line FHC exhibits tumorigenic phenotype, complex karyotype, and TP53 gene mutation. [Cancer Genet Cytogenet. 197: 107-16.](#)
3. Kalinina, T. *et al.* (2010) Establishment and characterization of a new human pancreatic adenocarcinoma cell line with high metastatic potential to the lung. [BMC Cancer.10: 295.](#)
4. Dallas, M.R. *et al.* (2012) Divergent roles of CD44 and carcinoembryonic antigen in colon cancer metastasis. [FASEB J. 226: 2648-56.](#)
5. Stern-Ginossar, N. *et al.* (2007) Intercellular Transfer of Carcinoembryonic Antigen from Tumor Cells to NK Cells. [J Immunol. 2007 Oct 1;179: 4424-34.](#)
6. Ferro, F. *et al.* (2011) Adipose tissue-derived stem cell in vitro differentiation in a three-dimensional dental bud structure. [Am J Pathol.178: 2299-310.](#)
7. Chao, A. *et al.* (2006) Molecular characterization of adenocarcinoma and squamous carcinoma of the uterine cervix using microarray analysis of gene expression. [Int J Cancer. 119: 91-8.](#)
8. Domenis, R. *et al.* (2015) Adipose tissue derived stem cells: in vitro and in vivo analysis of a standard and three commercially available cell-assisted lipotransfer techniques. [Stem Cell Res Ther. 6: 2.](#)
9. Wicklein, D. *et al.* (2018) CEACAM1 promotes melanoma metastasis and is involved in the regulation of the EMT associated gene network in melanoma cells. [Sci Rep. 8 \(1\): 11893.](#)
10. Caponnetto, F. *et al.* (2020) Human Adipose-Derived Stem Cells in Madelung's

---

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA1744F>  
10041

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA928F\)](#)

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M388022:210804'

Printed on 13 Feb 2024