

## Datasheet: MCA1750

**BATCH NUMBER 155871**

<b>Description:</b>	MOUSE ANTI PIG CD45RC
<b>Specificity:</b>	CD45RC
<b>Format:</b>	S/N
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	MIL15
<b>Isotype:</b>	IgM
<b>Quantity:</b>	2 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Pig
<b>Product Form</b>	Tissue Culture Supernatant - liquid
<b>Preparation</b>	Tissue Culture Supernatant containing 0.2M Tris/HCl pH7.4 and 5-10% foetal calf serum
<b>Preservative Stabilisers</b>	0.09% Sodium Azide
<b>Immunogen</b>	Cells isolated from porcine mesenteric lymph node.
<b>RRID</b>	AB_2174250

<b>Specificity</b>	<b>Mouse anti Pig CD45RC, clone MIL5</b> recognizes a restricted CD45 epitope encoded by the CD45 C exon and will bind and immunoprecipitate both the ~210 kDa RC and ~226 kDa RAC isoforms. MIL15, binds to an epitope alternative to that bound by clone <a href="#">MIL5</a> , also against CD45RC differing in their ability to stain B cells in tissue sections ( <a href="#">Zuckerman et al. 2001</a> ).
<b>Flow Cytometry</b>	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells in 100ul.
<b>References</b>	1. Zuckermann, F.A. <i>et al.</i> (2001) Characterization of monoclonal antibodies assigned to the CD45 subgroup of the Third International Swine CD Workshop. <a href="#">Vet Immunol Immunopathol. 80 (1-2): 165-74.</a>
<b>Further Reading</b>	1. Piriou-Guzylack, L. (2008) Membrane markers of the immune cells in swine: an update. <a href="#">Vet Res. 39: 54.</a>
<b>Storage</b>	Store at +4°C for one month or at -20°C for longer.  This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10053 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA1750">https://www.bio-rad-antibodies.com/SDS/MCA1750</a> 10053
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgM (STAR138...) [Alk. Phos.](#)

Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgM NEGATIVE CONTROL \(MCA692\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M365675:200529'

Printed on 29 Feb 2024