

Datasheet: MCA1750

BATCH NUMBER L1512

Description:	MOUSE ANTI PIG CD45RC
Specificity:	CD45RC
Format:	S/N
Product Type:	Monoclonal Antibody
Clone:	MIL15
Isotype:	IgM
Quantity:	2 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Tissue Culture Supernatant - liquid
Preparation	Tissue Culture Supernatant containing 0.2M Tris/HCl pH7.4 and 5-10% foetal calf serum
Preservative Stabilisers	0.09% Sodium Azide
Immunogen	Cells isolated from porcine mesenteric lymph node.
RRID	AB_2174250

Specificity	Mouse anti Pig CD45RC, clone MIL5 recognizes a restricted CD45 epitope encoded by the CD45 C exon and will bind and immunoprecipitate both the ~210 kDa RC and ~226 kDa RAC isoforms. MIL15, binds to an epitope alternative to that bound by clone MIL5 , also against CD45RC differing in their ability to stain B cells in tissue sections (Zuckerman et al. 2001).
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	1. Zuckermann, F.A. <i>et al.</i> (2001) Characterization of monoclonal antibodies assigned to the CD45 subgroup of the Third International Swine CD Workshop. Vet Immunol Immunopathol. 80 (1-2): 165-74.
Further Reading	1. Piriou-Guzylack, L. (2008) Membrane markers of the immune cells in swine: an update. Vet Res. 39: 54.
Storage	Store at +4°C for one month or at -20°C for longer. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10053 available at: https://www.bio-rad-antibodies.com/SDS/MCA1750 10053
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgM (STAR138...) [Alk. Phos.](#)

Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)

Recommended Negative Controls

[MOUSE IgM NEGATIVE CONTROL \(MCA692\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M365675:200529'

Printed on 29 Feb 2024