

Datasheet: MCA1773

Description:	MOUSE ANTI RAT CD61
Specificity:	CD61
Other names:	INTEGRIN BETA 3 CHAIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	F11
Isotype:	IgG1
Quantity:	0.25 mg

Product Details

RRID AB_322753

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Rat

Species Cross Reactivity Reacts weakly with:Human
N.B. Antibody reactivity and working conditions may vary between species.

Product Form Purified IgG - liquid

Buffer Solution Phosphate buffered saline

Preservative Stabilisers 0.09% Sodium Azide

Carrier Free Yes

Approx. Protein Concentrations IgG concentration 1.0 mg/ml

Immunogen	Bone cell suspension from SD rats.
Fusion Partners	Spleen cells from immunized Balb/c mice were fused with the cells of the mouse X63-Ag8 myeloma cell line.
Specificity	<p>Mouse anti Rat CD61 antibody, clone F11 recognizes the rat CD61 cell surface antigen, also known as the integrin b3 chain.</p> <p>CD61 is expressed by platelets, megakaryocytes and osteoclasts. Clone F11 has been shown to inhibit the <i>in vitro</i> attachment of osteoclasts to various proteins including vitronectin, fibrinogen and von Willebrand Factor (Sick et al. 2009).</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> Helfrich, M.H. <i>et al.</i> (1992) Integrins on rat osteoclasts: characterization of two monoclonal antibodies (F4 and F11) to rat beta 3. J Bone Miner Res. 7 (3): 345-51. Shankar, G. <i>et al.</i> (1993) Integrin receptor-mediated mobilisation of intranuclear calcium in rat osteoclasts. J Cell Sci. 105 (Pt 1): 61-8. Sick, E. <i>et al.</i> (2009) Activation of CD47 receptors causes histamine secretion from mast cells. Cell Mol Life Sci. 66: 1271-82. Wang, Z. <i>et al.</i> (2015) Microparticle levels after arterial injury and NO therapy in diabetes. J Surg Res. pii: S0022-4804(15)00867-7
Storage	<p>Store at +4°C for one month or at -20°C for longer.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL (MCA1209)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M341007:190109'

Printed on 11 Oct 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)