

Datasheet: MCA1782EL

Description:	MOUSE ANTI BOVINE INTERLEUKIN-12/23:Low Endotoxin
Specificity:	IL-12 / IL-23
Format:	Low Endotoxin
Product Type:	Monoclonal Antibody
Clone:	CC301
Isotype:	IgG2a
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	▪			1/10
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			8ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			1/50 - 1/100
Functional Assays	▪			1/100
ELISpot	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Membrane permeabilization is required for this application. The use of Leucoperm (Product Code [BUF09](#)) is recommended for this purpose.

Target Species

Bovine

Species Cross Reactivity

Reacts with: Human, Sheep, African Buffalo
Does not react with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	<0.01 EU/μg
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant bovine IL-12.
External Database Links	<p>UniProt:</p> <p>P54349 Related reagents</p> <p>P46282 Related reagents</p> <p>Entrez Gene:</p> <p>281856 IL12A Related reagents</p> <p>281857 IL12B Related reagents</p>
RRID	AB_616909
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Specificity	<p>Mouse anti Bovine Interleukin-12/23 antibody, clone CC301 recognizes the p40 subunit of bovine interleukin-12, and binds to both the free subunit and the intact heterodimer. The p40 subunit is also known as IL-12B and can form a heterodimer with either IL-12A or IL-23A.</p> <p>Mouse anti Bovine Interleukin-12/23 antibody, clone CC301 may be used as a capture antibody in an ELISA in combination with biotinylated Mouse anti Bovine Interleukin-12/23 antibody, clone CC326 (MCA2173B) as a detection reagent.</p>
Flow Cytometry	Use 10μl of the suggested working dilution to label 10 ⁶ cells in 100μl
ELISA	This antibody may be used as a capture antibody in a sandwich ELISA in combination with biotinylated clone CC326 (MCA2173B) as detection reagent (Bannerman et al. 2004).
References	<ol style="list-style-type: none"> Hope, J.C. <i>et al.</i> (2002) Development of detection methods for ruminant interleukin (IL)-12. J Immunol Methods. 266 (1-2): 117-26. Wenz, J.R. <i>et al.</i> (2010) Factors associated with concentrations of select cytokine and acute phase proteins in dairy cows with naturally occurring clinical mastitis. J Dairy Sci.

[93: 2458-70.](#)

3. Rinaldi, M. *et al.* (2010) A sentinel function for teat tissues in dairy cows: dominant innate immune response elements define early response to *E. coli* mastitis. [Funct Integr Genomics. 10: 21-38.](#)
4. Ferret-Bernard, S. *et al.* (2011) Mesenteric lymph node cells from neonates present a prominent IL-12 response to CpG oligodeoxynucleotide via an IL-15 feedback loop of amplification. [Vet Res. 42:19.](#)
5. Contreras, V. *et al.* (2010) Existence of CD8 α -like dendritic cells with a conserved functional specialization and a common molecular signature in distant mammalian species. [J Immunol. 185: 3313-25.](#)
6. Bannerman, D.D. *et al.* (2004) *Escherichia coli* and *Staphylococcus aureus* elicit differential innate immune responses following intramammary infection. [Clin Diagn Lab Immunol. 11: 463-72.](#)
7. Davis, T.L. and Pate, J.L. (2007) Bovine luteal cells stimulate proliferation of major histocompatibility nonrestricted gamma delta T cells. [Biol Reprod. 77: 914-22.](#)
8. Ferret-Bernard, S. (2010) Cellular and molecular mechanisms underlying the strong neonatal IL-12 response of lamb mesenteric lymph node cells to R-848. [PLoS One. 5: e13705.](#)
9. Shoda, L.K. *et al.* (2001) DNA from protozoan parasites *Babesia bovis*, *Trypanosoma cruzi*, and *T. brucei* is mitogenic for B lymphocytes and stimulates macrophage expression of interleukin-12, tumor necrosis factor alpha, and nitric oxide. [Infect Immun. 69: 2162-71.](#)
10. Souza, M. *et al.* (2008) Pathogenesis and immune responses in gnotobiotic calves after infection with the genogroup II.4-HS66 strain of human norovirus. [J Virol. 82: 1777-86.](#)
11. Verhelst, D. *et al.* (2014) Parasite distribution and associated immune response during the acute phase of *Toxoplasma gondii* infection in sheep. [BMC Vet Res. 10: 293.](#)
12. Beechler BR *et al.* (2015) Enemies and turncoats: bovine tuberculosis exposes pathogenic potential of Rift Valley fever virus in a common host, African buffalo (*Syncerus caffer*). [Proc Biol Sci. 282 \(1805\): .](#)
13. Pomeroy B *et al.* (2015) Monocyte-derived dendritic cells from late gestation cows have an impaired ability to mature in response to *E. coli* stimulation in a receptor and cytokine-mediated fashion. [Vet Immunol Immunopathol. 167 \(1-2\): 22-9.](#)
14. Stabel, J.R. & Bannantine, J.P. (2021) Reduced tissue colonization of Mycobacterium avium subsp. paratuberculosis in neonatal calves vaccinated with a cocktail of recombinant proteins. [Vaccine. 39 \(23\): 3131-40.](#)
15. Rodrigues, V. *et al.* (2017) Development of a bead-based multiplexed assay for simultaneous quantification of five bovine cytokines by flow cytometry. [Cytometry A. 91 \(9\): 901-7.](#)
16. Stabel, J.R. *et al.* (2021) Comparative cellular immune responses in calves after infection with Mycobacterium avium subsp. paratuberculosis, M. avium subsp. avium, M. kansasii and M. bovis. [Vet Immunol Immunopathol. 237: 110268.](#)
17. Ciliberti, M.G. *et al.* (2020) Nexus Between Immune Responses and Oxidative Stress: The Role of Dietary Hydrolyzed Lignin in *ex vivo* Bovine Peripheral Blood Mononuclear Cell Response. [Front Vet Sci. 7: 9.](#)
18. Tavalire, H.F. *et al.* (2019) Risk alleles for tuberculosis infection associate with reduced immune reactivity in a wild mammalian host. [Proc Biol Sci. 286 \(1907\): 20190914.](#)
19. Stabel, J.R. *et al.* (2020) Comparison of Sheep, Goats, and Calves as Infection

Models for *Mycobacterium avium* subsp. *paratuberculosis*. [Vet Immunol Immunopathol. 225: 110060.](#)

20. Ihedioha, O. *et al.* (2020) Poor stimulation of bovine dendritic cells by *Mycobacterium bovis* culture supernatant and surface extract is associated with decreased activation of ERK and NF- κ B and higher expression of SOCS1 and 3. [Innate Immun. 26 \(6\): 537-546.](#)

Storage	Store at -20°C only. This product should be stored undiluted. Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/MCA1782EL 10162
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL:Low Endotoxin \(MCA929EL\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M411493:221103'

Printed on 26 Feb 2024