

## Datasheet: MCA1815T

**BATCH NUMBER 162025**

<b>Description:</b>	MOUSE ANTI HUMAN CD68
<b>Specificity:</b>	CD68
<b>Other names:</b>	MACROSIALIN
<b>Format:</b>	S/N
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	514H12
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin (1)	▪			1/40 - 1/80
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting		▪		
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Tris/EDTA buffer pH 9.0 is recommended for this purpose. *N.B. The epitope recognised by clone 514H12 is sensitive to peroxide. It is therefore recommended that any peroxidase-blocking steps be performed after incubation with this antibody.***

<b>Target Species</b>	Human
<b>Product Form</b>	Tissue Culture Supernatant - liquid

<b>Preservative Stabilisers</b>	<0.1% Sodium Azide (NaN <sub>3</sub> )
<b>Immunogen</b>	Fusion protein corresponding to external domain of human CD68.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P34810</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">968</a>    CD68    <a href="#">Related reagents</a></p>
<b>RRID</b>	AB_2074721
<b>Fusion Partners</b>	Spleen cells from immunized mice were fused with cells of the mouse p3-NS1-Ag4-1 myeloma cell line
<b>Specificity</b>	<b>Mouse anti Human CD68 antibody, clone 514H12</b> recognizes the human CD68 cell surface antigen, a ~110 kDa glycoprotein primarily expressed by macrophages and monocytes.
<b>Histology Positive Control Tissue</b>	Human Tonsil
<b>References</b>	<ol style="list-style-type: none"> <li>Hameed, A. <i>et al.</i> (1994) Immunohistochemical expression of CD68 antigen in human peripheral blood T cells. <a href="#">Hum Pathol. 25 (9): 872-6.</a></li> <li>Madden, L.R. <i>et al.</i> (2010) Proangiogenic scaffolds as functional templates for cardiac tissue engineering. <a href="#">Proc Natl Acad Sci U S A. 107 (34): 15211-6.</a></li> <li>Muthana, M. <i>et al.</i> (2011) Use of macrophages to target therapeutic adenovirus to human prostate tumors. <a href="#">Cancer Res. 71 (5): 1805-15.</a></li> <li>da Costa, C.E. <i>et al.</i> (2005) Presence of osteoclast-like multinucleated giant cells in the bone and nonostotic lesions of Langerhans cell histiocytosis. <a href="#">J Exp Med. 201 (5): 687-93.</a></li> <li>Wang, X. <i>et al.</i> (2006) Monocyte/macrophage and T-cell infiltrates in peritoneum of patients with ovarian cancer or benign pelvic disease. <a href="#">J Transl Med. 4: 30.</a></li> <li>Coury, F. <i>et al.</i> (2008) Langerhans cell histiocytosis reveals a new IL-17A-dependent pathway of dendritic cell fusion. <a href="#">Nat Med. 14: 81-7.</a></li> <li>Rodriguez-Agudo, D. <i>et al.</i> (2006) Localization of StarD5 cholesterol binding protein. <a href="#">J Lipid Res. 47: 1168-75.</a></li> <li>Achard, V. <i>et al.</i> (2007) Renin receptor expression in human adipose tissue. <a href="#">Am J Physiol Regul Integr Comp Physiol. 292: R274-82.</a></li> <li>Paulmyer-Lacroix, O. <i>et al.</i> (2006) Expression of adrenomedullin in adipose tissue of lean and obese women. <a href="#">Eur J Endocrinol. 155: 177-85.</a></li> <li>Hever, A. <i>et al.</i> (2007) Human endometriosis is associated with plasma cells and overexpression of B lymphocyte stimulator. <a href="#">Proc Natl Acad Sci U S A. 104: 12451-6.</a></li> <li>Silaghi, A. <i>et al.</i> (2007) Expression of adrenomedullin in human epicardial adipose tissue: role of coronary status. <a href="#">Am J Physiol Endocrinol Metab. 293: E1443-50.</a></li> <li>Kaibara, N. <i>et al.</i> (2008) Comparative histopathological analysis between tenosynovitis and joint synovitis in rheumatoid arthritis. <a href="#">Histopathology. 52: 856-64.</a></li> <li>Moskovszky, L. <i>et al.</i> (2009) Genomic instability in giant cell tumor of bone. A study of</li> </ol>

52 cases using DNA ploidy, relocalization FISH, and array-CGH analysis. [Genes Chromosomes Cancer. 48: 468-79.](#)

14. Angel, C.E. *et al.* (2007) CD14+ antigen-presenting cells in human dermis are less mature than their CD1a+ counterparts. [Int Immunol. 19: 1271-9.](#)

15. de Vos van Steenwijk PJ *et al.* (2013) Tumor-infiltrating CD14-positive myeloid cells and CD8-positive T-cells prolong survival in patients with cervical carcinoma. [Int J Cancer. 133 \(12\): 2884-94.](#)

16. Quispel, W.T. *et al.* (2016) Tertiary lymphoid structures are confined to patients presenting with unifocal Langerhans Cell Histiocytosis. [Oncoimmunology. 5 \(8\): e1164364.](#)

17. Toh, M.L. *et al.* (2014) Bone- and cartilage-protective effects of a monoclonal antibody against colony-stimulating factor 1 receptor in experimental arthritis. [Arthritis Rheumatol. 66 \(11\): 2989-3000.](#)

18. Schirmer, L. *et al.* (2014) Differential loss of KIR4.1 immunoreactivity in multiple sclerosis lesions. [Ann Neurol. 75 \(6\): 810-28.](#)

19. Schrevel, M. *et al.* (2017) Autocrine expression of the epidermal growth factor receptor ligand heparin-binding EGF-like growth factor in cervical cancer. [Int J Oncol. 50 \(6\): 1947-54.](#)

20. Zhang, W. *et al.* (2013) Myeloid clusters are associated with a pro-metastatic environment and poor prognosis in smoking-related early stage non-small cell lung cancer. [PLoS One. 8 \(5\): e65121.](#)

21. Bobosha, K. *et al.* (2014) T-cell regulation in lepromatous leprosy. [PLoS Negl Trop Dis. 8 \(4\): e2773.](#)

22. van Poppel, P.C.M. *et al.* (2018) Inflammation in the subcutaneous adipose tissue does not attenuate endothelial function in subjects with diabetes mellitus and subjects with dyslipidaemia and hypertension: A cross-sectional study. [Endocrinol Diabetes Metab. 1 \(3\): e00020.](#)

23. Yuan, A. *et al.* (2020) Histologic analysis of medication-related osteonecrosis of the jaw compared with antiresorptive-exposed bone and other infectious, inflammatory, and necrotic jaw diseases. [Oral Surg Oral Med Oral Pathol Oral Radiol. 129 \(2\): 133-40.](#)

24. Durán, D. *et al.* (2020) Histological and Immunohistochemical Study of Wounds in Sheep Skin in Maggot Therapy by Using *Protophormia terraenovae*. (Diptera: Calliphoridae.) Larvae. [J Med Entomol. 57 \(2\): 369-76.](#)

---

<b>Storage</b>	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
----------------	---

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10055 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA1815T">https://www.bio-rad-antibodies.com/SDS/MCA1815T</a> 10055
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> , <a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M389336:210806'

Printed on 29 Feb 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)