

Datasheet: MCA1818

BATCH NUMBER 173102

Description:	MOUSE ANTI HUMAN MCP-3
Specificity:	MCP-3
Other names:	CCL7
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	h.mcp.3
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/1000
Immunoprecipitation			▪	
Western Blotting			▪	
Functional Assays	▪			1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilized
Reconstitution	Reconstitute with 0.5 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently vortexed after reconstitution and microcentrifuged before use.
Preparation	Purified IgG prepared by affinity chromatography on Protein A from ascites

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant human MCP-3 (Macrophage Chemoattractant Protein-3).
External Database Links	<p>UniProt: P80098 Related reagents</p> <p>Entrez Gene: 6354 CCL7 Related reagents</p>
Synonyms	MCP3, SCYA6, SCYA7
RRID	AB_2259751
Fusion Partners	Spleen cells from immunised C3H/HeJ mice were fused with cells of a mouse myeloma cell line.
Specificity	Mouse anti Human MCP-3 antibody, clone h.mcp.3 recognizes human C-C motif chemokine 7, also known as Monocyte chemoattractant protein 3, MCP-3, CCL7, NC28 or Small-inducible cytokine A7. MCP-3 is a 99 amino acid (including signal peptide sequence) ~9kDa chemotactic glycoprotein for monocytes and eosinophils. No cross reactivity is observed with MCP-1, MCP-2 or other structurally related chemokines.
References	1. Al-Bannawi, A. <i>et al.</i> (2011) <i>Chlamydia pneumoniae</i> induces chemokine expression by platelets in patients with atherosclerosis. Med Princ Pract. 20 (5): 438-43.
Further Reading	1. Opendakker, G. <i>et al.</i> (1993) Human monocyte chemotactic protein-3 (MCP-3): molecular cloning of the cDNA and comparison with other chemokines. Biochem Biophys Res Commun. 191: 535-42.
Storage	This product is shipped at ambient temperature. Prior to reconstitution store at -20°C. After reconstitution store at -20°C. This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety	Material Safety Datasheet documentation #10294 available at:

Information <https://www.bio-rad-antibodies.com/SDS/MCA1818>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M440258:250523'

Printed on 19 Mar 2026