

Datasheet: MCA1858B

BATCH NUMBER 155517

Description:	MOUSE ANTI M13 BACTERIOPHAGE:Biotin
Specificity:	M13 BACTERIOPHAGE
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	E1
Isotype:	IgG3
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1ug/ml - 5ug/ml
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Viral
Product Form	Purified IgG conjugated to Biotin - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.08% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	M13 Phage coat proteins.

RRID	AB_323105
Specificity	<p>Mouse anti M13 Bacteriophage antibody, clone E1 recognizes the gp3 (G3P) phage coat protein of the M13 bacteriophage widely used in the preparation of 'phage-display' immunoglobulin libraries. G3P plays a role in the entry of the viral genome into the bacterial host..</p> <p>This antibody may be used to detect binding of ScFv antibodies developed in the M13 system.</p>
References	1. Jaye, D.L. <i>et al.</i> (2003) Use of real-time polymerase chain reaction to identify cell- and tissue-type-selective peptides by phage display. Am J Pathol. 162:1419-29.
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA1858B 10040
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M350396:190307'

Printed on 18 Jan 2024