

Datasheet: MCA1858R

BATCH NUMBER 149753

Description:	MOUSE ANTI M13 BACTERIOPHAGE
Specificity:	M13 BACTERIOPHAGE
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	E1
Isotype:	IgG3
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			10ug/ml
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Viral
Product Form	Purified IgG - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	0.5 mg/ml
Immunogen	M13 Phage coat proteins.

Specificity **Mouse anti M13 Bacteriophage antibody, clone E1** recognizes the gp3 (G3P) phage coat protein of the M13 bacteriophage widely used in the preparation of 'phage-display' immunoglobulin libraries. G3P plays a role in the entry of the viral genome into the bacterial host..

This antibody may be used to detect binding of ScFv antibodies developed in the M13 system.

Flow Cytometry Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul

References 1. Jaye, D.L.*et al.* (2003) Use of real-time polymerase chain reaction to identify cell- and tissue-type-selective peptides by phage display. [Am J Pathol. 162:1419-29.](#)

Storage Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA1858R>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)

Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M365830:200529'

Printed on 18 Jan 2024