

Datasheet: MCA1876

Description:	MOUSE ANTI HUMAN CD147
Specificity:	CD147
Other names:	NEUROTHELIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	MEM-M6/1
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting (1)	▪			

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)MEM-M6/1 recognizes CD147 under non-reducing conditions.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by Protein A affinity chromatography.
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1 mg/ml
Immunogen	CD147Rg, consisting of the cDNA coding for the extracellular region of human CD147.
External Database Links	<p>UniProt: P35613 Related reagents</p> <p>Entrez Gene: 682 BSG Related reagents</p>
RRID	AB_322225
Fusion Partners	Spleen cells from immunized Balb/c mice were fused with SP2/0 mouse myeloma cells.
Specificity	<p>Mouse anti Human CD147 antibody, clone MEM-M6/1 recognizes the human CD147 cell surface antigen, also known as Basigin, EMMPRIN and collagenase stimulatory factor. CD147 is a 385 amino acid single pass type 1 trans-membrane glycoprotein bearing 3 potential N-glycosylation sites, a single N-terminal (distal) Ig-like C2-type domain and a more proximal Ig-like V-type domain in its extracellular region. Mouse anti Human CD147 antibody, clone MEM-M6/1 was raised and screened against recombinant extracellular region of huCD147 expressed in CHO cells, along with a number of other MEM-M6 clones. Epitope mapping indicates that clone MEM-M6/1 binds to an epitope within the distal Ig-like C-2 type domain. This is confirmed by flow cytometry using beads coated with constructs containing only the distal Ig-like or proximal Ig-like domains (Koch et al. 1999).</p> <p>Mouse anti Human CD147 antibody, clone MEM-M6/1 binds to both un-stimulated and phytohemagglutinin (PHA) stimulated (activated) T lymphocytes with levels of CD147 enhanced following PHA stimulation. CD147 is also expressed by peripheral blood and endothelial cells as well as many cultured cells of hematopoietic and non-hematopoietic origin (Koch et al. 1999). In addition to clone MEM-M6/1 we are also able to offer the MEM-M6/6 clone derived from the same fusion which recognizes an epitope in the membrane proximal Ig-like V type domain and has inhibitory effects on CD3 induced T cell activation (Kock et al. 1999)</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells or 100ul whole blood.
References	<ol style="list-style-type: none"> 1. Koch, C. <i>et al.</i> (1999) T cell activation-associated epitopes of CD147 in regulation of the T cell response, and their definition by antibody affinity and antigen density. Int Immunol. 11 (5): 777-86. 2. Yang, Y. <i>et al.</i> (2008) Cyclophilin A up-regulates MMP-9 expression and adhesion of monocytes/macrophages via CD147 signalling pathway in rheumatoid arthritis. Rheumatology (Oxford). 47: 1299-310. 3. Schneiderhan, W. <i>et al.</i> (2007) Pancreatic stellate cells are an important source of

- MMP-2 in human pancreatic cancer and accelerate tumor progression in a murine xenograft model and CAM assay. [J Cell Sci. 120: 512-9.](#)
4. Chen, Y. *et al.* (2009) Upregulation of HAb18G/CD147 in activated human umbilical vein endothelial cells enhances the angiogenesis. [Cancer Lett. 278: 113-21.](#)
 5. Schneiderhan, W. *et al.* (2009) CD147 silencing inhibits lactate transport and reduces malignant potential of pancreatic cancer cells in in vivo and in vitro models [Gut. 58: 1391-8.](#)
 6. Wright, G.J. *et al.* (2014) Immunogenic compositions and expression systems. [Patent Publication number. US20140093540 A1](#)
 7. Huang, Z. *et al.* (2013) Overexpression of CD147 contributes to the chemoresistance of head and neck squamous cell carcinoma cells. [J Oral Pathol Med. 42: 541-6.](#)
 8. Hu, J. *et al.* (2010) Involvement of HAb18G/CD147 in T cell activation and immunological synapse formation. [J Cell Mol Med. 14 \(8\): 2132-43.](#)
 9. Bernard, S.C. *et al.* (2014) Pathogenic *Neisseria meningitidis* utilizes CD147 for vascular colonization. [Nat Med. 20 \(7\): 725-31.](#)
 10. Nassif, X. *et al.* (2015) Cd147 as receptor for pilus-mediated adhesion of meningococci to vascular endothelia. [Patent Publication number US20150110806 A1](#)
 11. Trakarnsanga, K. *et al.* (2017) An immortalized adult human erythroid line facilitates sustainable and scalable generation of functional red cells. [Nat Commun. 8: 14750.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements PrecisionAb is a trademark of Bio-Rad Laboratories

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

- | | |
|---|--|
| Goat Anti Mouse IgG (STAR77...) | HRP |
| Rabbit Anti Mouse IgG (STAR12...) | RPE |
| Rabbit Anti Mouse IgG (STAR8...) | DyLight@800 |
| Goat Anti Mouse IgG (STAR70...) | FITC |
| Goat Anti Mouse IgG IgA IgM (STAR87...) | Alk. Phos. , HRP |
| Goat Anti Mouse IgG (Fc) (STAR120...) | FITC , HRP |
| Rabbit Anti Mouse IgG (STAR9...) | FITC |

Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

Goat Anti Mouse IgG (STAR76...) [RPE](#)

Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M393633:220113'

Printed on 23 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)