

Datasheet: MCA1893

BATCH NUMBER 165721

Description:	MOUSE ANTI CANINE DISTEMPER VIRUS
Specificity:	CANINE DISTEMPER VIRUS
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	DV2-12
Isotype:	IgG2b
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin (1)	▪			
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Serotec target unmasking fluid (STUF) product code [BUF025A/B/C](#) is recommended for this purpose.

Target Species	Viral
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
RRID	AB_323083
Specificity	<p>Mouse anti canine distemper virus antibody, clone DV2-12 recognizes canine distemper virus nucleoprotein. In western blots the antibody recognizes bands of ~76 kDa, 68 kDa, 54 kDa and 32 kDa under reducing conditions.</p> <p>Mouse anti canine distemper virus antibody, clone DV2-12 is useful for the detection of CDV in immunofluorescence procedures, in which the staining seen has a cytoplasmic pattern in infected cells. It is also effective in paraffin embedded material.</p> <p>Mouse anti canine distemper virus antibody, clone DV2-12 is anti morbillivirus specific and reacts with phocine distemper virus, raccoon distemper virus and human measles virus. It does not cross react with CAV2, CCV, CPI or CPV.</p>
References	<ol style="list-style-type: none"> Liang, C.T. <i>et al.</i> (2007) A non-biotin polymerized horseradish-peroxidase method for the immunohistochemical diagnosis of canine distemper. J Comp Pathol. 136 (1): 57-64. Sánchez-Migallón Guzmán, D. <i>et al.</i> (2008) Aleutian disease serology, protein electrophoresis, and pathology of the European mink (<i>Mustela lutreola</i>) from Navarra, Spain. J Zoo Wildl Med. 39: 305-13. Bathen-Noethen, A. <i>et al.</i> (2008) Magnetic resonance imaging findings in acute canine distemper virus infection. J Small Anim Pract. 49: 460-7. Guvenc, T. <i>et al.</i> (2008) Immunohistochemical distribution of alpha B-crystallin in the cerebellum of dogs infected with canine distemper virus. Acta Vet Hung. 56: 117-23. Headley, S.A. and Sukura, A. (2009) Naturally occurring systemic canine distemper virus infection in a pup Braz J Vet Pathol 2: 95-101. Headley, S.A. <i>et al.</i> (2009) Molecular detection of Canine distemper virus and the immunohistochemical characterization of the neurologic lesions in naturally occurring old dog encephalitis. J Vet Diagn Invest. 21: 588-97. Headley, S.A. <i>et al.</i> (2009) Diagnostic exercise: Tyzzer's disease, distemper, and coccidiosis in a pup. Vet Pathol. 46: 151-4. Spitzbarth, I. <i>et al.</i> (2010) Immunohistochemical characterization of inflammatory and glial responses in a case of necrotizing leucoencephalitis in a French bulldog. J Comp Pathol. 142 (2-3): 235-41. Woo, G.H. <i>et al.</i> (2010) Canine distemper virus infection in fennec fox (<i>Vulpes zerda</i>). J Vet Med Sci. 72 (8): 1075-9. Kubo, M. <i>et al.</i> (2010) Meningoencephalitis associated with Sarcocystis spp. in a free-living Japanese raccoon dog (<i>Nyctereutes procyonoides viverrinus</i>). J Comp Pathol. 143: 185-9. Kyöstilä, K. <i>et al.</i> (2012) A SEL1L mutation links a canine progressive early-onset cerebellar ataxia to the endoplasmic reticulum-associated protein degradation (ERAD) machinery. PLoS Genet. 8 (6): e1002759. Szentiks, C.A. <i>et al.</i> (2013) Polar bear encephalitis: establishment of a comprehensive

next-generation pathogen analysis pipeline for captive and free-living wildlife. [J Comp Pathol. 2014 May;150\(4\):474-88.](#)

13. Rentería-Solís, Z. *et al.* (2014) Canine distemper outbreak in raccoons suggests pathogen interspecies transmission amongst alien and native carnivores in urban areas from Germany. [Vet Microbiol. 174 \(1-2\): 50-9.](#)

14. Bak, E.J. & , Woo, G.H. (2016) Canine Distemper in Albino Ferrets (*Mustela putorius furo*). [J. Comp Path. 154 \(1\): 121.](#)

15. Çomaklı, S. *et al.* (2020) Canine distemper virus induces downregulation of GABA_A, GABA_B, and GAT1 expression in brain tissue of dogs. [Arch Virol. 165 \(6\): 1321-31.](#)

16. Iribarnegaray, V. *et al.* (2024) Droplet Digital PCR Enhances Sensitivity and Accuracy in Canine Distemper Virus Detection [Preprints 18 Sept \[Epub ahead of print\].](#)

17. Garcia, P.A.T. *et al.* (2022) Distribution of canine distemper virus and nectin-4 in raccoon (*Procyon lotor*) skin. [Vet Pathol. 59 \(5\): 782-786.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA1893>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),
[FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M419062:230516'

Printed on 20 Sep 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)