

Datasheet: MCA1896B

Description:	MOUSE ANTI DOG IgE:Biotin
Specificity:	IgE
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	E6-71A1
Isotype:	IgG1
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Dog
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Affinity purified IgE preparation from heavily parasitized and allergic dog serum.

RRID AB_323085

Specificity **Mouse anti Dog IgE antibody, clone E6-71A1** recognises canine IgE. E6-2A1 does not cross react with dog IgM, IgA, IgG1 or IgG2. Western blot analysis against affinity purified dog IgE using Mouse anti Dog IgE clone E6-71A1 demonstrates a single major band of ~62-65 kDa under reducing conditions.

ELISA The antibody may be used as a detection antibody in a sandwich ELISA in combination with MCA1895 as the capture antibody.

Western Blotting Under reducing conditions MCA1896B detects a band of approximately 62-65 kDa in samples of affinity purified dog IgE.

References

1. Ognjenovic, J. *et al.* (2013) Immunoproteomic characterization of *Ambrosia artemisiifolia* pollen allergens in canine atopic dermatitis. [Vet Immunol Immunopathol. 155: 38-47.](#)
2. Martins, L.M. *et al.* (2017) Allergy to grass pollen: mapping of *Dactylis glomerata* and *Phleum pratense* allergens for dogs by two-dimensional immunoblotting. [Postepy Dermatol Alergol. 34 \(1\): 60-9.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M389349:210806'

Printed on 21 Mar 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)