

## Datasheet: MCA1971SBV570

<b>Description:</b>	MOUSE ANTI PIG CD16:StarBright Violet 570
<b>Specificity:</b>	CD16
<b>Other names:</b>	FcRIII
<b>Format:</b>	StarBright Violet 570
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	G7
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 TESTS/0.5ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Pig		
<b>Product Form</b>	Purified IgG conjugated to StarBright Violet 570 - liquid		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	StarBright Violet 570	404	571
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ) 1% Bovine Serum Albumin 0.1% Pluronic F68 0.1% PEG 3350 0.05% Tween 20		

<b>Approx. Protein Concentrations</b>	For information on the concentration of our StarBright Dye conjugated reagents please visit our <a href="#">FAQ</a> page.
<b>Immunogen</b>	Porcine peripheral blood leucocytes
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">Q28942</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">397684</a>    FCGR3B    <a href="#">Related reagents</a></p>
<b>Fusion Partners</b>	Spleen cells from immunized Balb/c mice were fused with cells of the mouse P3-X63-Ag8.653 myeloma cell line
<b>Specificity</b>	<p><b>Mouse anti Pig CD16, clone G7</b> recognizes porcine CD16 also known as Fc-gamma RIII or the low affinity IgG (Fc) receptor III. Clone G7 was clustered as CD16 at the Second International Workshop to Define Swine Cluster of Differentiation (CD) Antigens (<a href="#">Saalmuller et al. 1998</a>).</p> <p>Mouse anti pig CD16 immunoprecipitates a protein of ~40 kDa from porcine neutrophils and NK cells (<a href="#">Wierda et al. 1993</a>). Subsequent cloning and characterization of the G7 molecule indicated that G7 was the porcine homologue of Human CD16 (<a href="#">Halloran et al. 1994</a>).</p>
<b>Flow Cytometry</b>	Use 5µl of the suggested working dilution to label 0.5x10 <sup>6</sup> cells in 100µl. Best practices suggest a 5 min centrifugation at 6,000g prior to sample application.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Wierda, W.G. <i>et al.</i> (1993) Two distinct porcine natural killer lytic trigger molecules as PNK-E/G7 molecular complex. <a href="#">Cell Immunol. 146 (2): 270-83.</a></li> <li>2. Halloran, P.J. <i>et al.</i> (1994) Biochemical characterization of the porcine Fc gamma RIII alpha homologue G7. <a href="#">Cell Immunol. 158 (2): 400-13.</a></li> <li>3. Sánchez, C. <i>et al.</i> (1999) The porcine 2A10 antigen is homologous to human CD163 and related to macrophage differentiation. <a href="#">J Immunol. 162 (9): 5230-7.</a></li> <li>4. Terzic, S. <i>et al.</i> (2002) Immunophenotyping of leukocyte subsets in peripheral blood and palatine tonsils of prefattening pigs. <a href="#">Vet Res Commun. 26: 273 - 83.</a></li> <li>5. Vincent, I.E. <i>et al.</i> (2003) Dendritic cells harbor infectious porcine circovirus type 2 in the absence of apparent cell modulation or replication of the virus. <a href="#">J Virol. 77: 13288 - 300.</a></li> <li>6. Summerfield, A. <i>et al.</i> (2003) Porcine peripheral blood dendritic cells and natural interferon-producing cells. <a href="#">Immunology 110: 440-9.</a></li> <li>7. Inman, C.F. <i>et al.</i> (2010) Rearing environment affects development of the immune system in neonates. <a href="#">Clin Exp Immunol. 160 (3): 431-9.</a></li> <li>8. Inman, C.F. <i>et al.</i> (2010) Dendritic cells interact with CD4 T cells in intestinal mucosa. <a href="#">J Leukoc Biol. 88 (3): 571-8.</a></li> <li>9. Devriendt, B. <i>et al.</i> (2010) Targeting of <i>Escherichia coli</i> F4 fimbriae to Fcgamma receptors enhances the maturation of porcine dendritic cells. <a href="#">Vet Immunol Immunopathol. 135 (3-4): 188-98.</a></li> </ol>

10. Gimeno, M. *et al.* (2011) Cytokine profiles and phenotype regulation of antigen presenting cells by genotype-I porcine reproductive and respiratory syndrome virus isolates. [Vet Res. 42: 9.](#)
11. Mussá, T. *et al.* (2011) Interaction of porcine conventional dendritic cells with swine influenza virus. [Virology 420: 125-34.](#)
12. Lecours, M.P. *et al.* (2011) Characterization of porcine dendritic cell response to *Streptococcus suis*. [Vet Res. 42: 72.](#)
13. Mair, K.H. *et al.* (2012) NKp46 expression discriminates porcine NK cells with different functional properties. [Eur J Immunol. 42: 1261-71.](#)
14. Hester, S.N. *et al.* (2012) Intestinal and systemic immune development and response to vaccination are unaffected by dietary (1,3/1,6)- $\beta$ -D-glucan supplementation in neonatal piglets. [Clin Vaccine Immunol. 19 \(9\): 1499-508.](#)
15. Thierry, A. *et al.* (2012) Identification of invariant natural killer T cells in porcine peripheral blood. [Vet Immunol Immunopathol. 149 \(3-4\): 272-9.](#)
16. Inman, C.F. *et al.* (2012) Neonatal colonisation expands a specific intestinal antigen-presenting cell subset prior to CD4 T-cell expansion, without altering T-cell repertoire. [PLoS One 7: e33707.](#)
17. Masure, D. *et al.* (2013) A Role for Eosinophils in the Intestinal Immunity against Infective *Ascaris suum* Larvae. [PLoS Negl Trop Dis. 7: e2138.](#)
18. Mair, K.H. *et al.* (2013) Porcine CD8 $\alpha$ dim<sup>-</sup>/NKp46<sup>high</sup> NK cells are in a highly activated state. [Vet Res. 44: 13.](#)
19. Kyrova, K. *et al.* (2014) The response of porcine monocyte derived macrophages and dendritic cells to *Salmonella typhimurium* and lipopolysaccharide. [BMC Vet Res. 10: 244.](#)
20. Waide, E.H. *et al.* (2015) Not All SCID Pigs Are Created Equally: Two Independent Mutations in the Artemis Gene Cause SCID in Pigs. [J Immunol. 195 \(7\): 3171-9.](#)
21. Auray, G. *et al.* (2016) Characterization and Transcriptomic Analysis of Porcine Blood Conventional and Plasmacytoid Dendritic Cells Reveals Striking Species-Specific Differences. [J Immunol. 197 \(12\): 4791-806.](#)
22. Suzuki, S. *et al.* (2016) Generation and characterization of RAG2 knockout pigs as animal model for severe combined immunodeficiency. [Vet Immunol Immunopathol. 178: 37-49.](#)
23. LeLuduec, J.B. *et al.* (2016) Intradermal vaccination with un-adjuvanted sub-unit vaccines triggers skin innate immunity and confers protective respiratory immunity in domestic swine. [Vaccine. 34 \(7\): 914-22.](#)
24. Waddell, L.A. *et al.* (2018) ADGRE1 (EMR1, F4/80) Is a Rapidly-Evolving Gene Expressed in Mammalian Monocyte-Macrophages. [Front Immunol. 9: 2246.](#)
25. Fernández-Caballero, T. *et al.* (2018) Phenotypic and functional characterization of porcine bone marrow monocyte subsets. [Dev Comp Immunol. 81: 95-104.](#)
26. Loss, H. *et al.* (2018) Effects of a pathogenic ETEC strain and a probiotic *Enterococcus faecium* strain on the inflammasome response in porcine dendritic cells. [Vet Immunol Immunopathol. 203: 78-87.](#)
27. Skovdal, S.M. *et al.* (2019) Inhaled nebulized glatiramer acetate against Gram-negative bacteria is not associated with adverse pulmonary reactions in healthy, young adult female pigs. [PLoS One. 14 \(10\): e0223647.](#)
28. Ferret-Bernard, S. *et al.* (2020) Maternal Supplementation of Food Ingredient (Prebiotic) or Food Contaminant (Mycotoxin) Influences Mucosal Immune System in Piglets. [Nutrients. 12 \(7\): 2115.](#)

29. Boettcher, A.N. *et al.* (2020) CD3ε(+) Cells in Pigs With Severe Combined Immunodeficiency Due to Defects in ARTEMIS. [Front Immunol. 11: 510.](#)
30. Van der Weken, H. *et al.* (2021) Antibody-Mediated Targeting of Antigens to Intestinal Aminopeptidase N Elicits Gut IgA Responses in Pigs. [Front Immunol. 12: 753371.](#)
31. Teuben, M. *et al.* (2021) Instant intra-operative neutropenia despite the emergence of banded (CD16<sup>dim</sup>/CD62L<sup>bright</sup>) neutrophils in peripheral blood - An observational study during extensive trauma-surgery in pigs. [Injury. 52 \(3\): 426-33.](#)
32. Teuben, M.P.J. *et al.* (2022) Standardized porcine unilateral femoral nailing is associated with changes in PMN activation status, rather than aberrant systemic PMN prevalence. [Eur J Trauma Emerg Surg. 48 \(3\): 1601-11.](#)
33. Zhao, H. *et al.* (2022) Development of *RAG2*<sup>-/-</sup> *IL2Rγ*<sup>-/-</sup> immune deficient FAH-knockout miniature pig. [Front Immunol. 13: 950194.](#)
34. Štěpánová, H. *et al.* (2022) Characterization of Porcine Monocyte-Derived Macrophages Cultured in Serum-Reduced Medium. [Biology \(Basel\). 11\(10\):1457.](#)
35. Monguíó-Tortajada, M. *et al.* (2022) Acellular cardiac scaffolds enriched with MSC-derived extracellular vesicles limit ventricular remodelling and exert local and systemic immunomodulation in a myocardial infarction porcine model. [Theranostics. 12 \(10\): 4656-70.](#)
36. Haach, V. *et al.* (2023) A polyvalent virosomal influenza vaccine induces broad cellular and humoral immunity in pigs. [Virology. 20 \(1\): 181.](#)
37. Álvarez, B. *et al.* (2023) Porcine Macrophage Markers and Populations: An Update. [Cells. 12 \(16\): 2103.](#)
38. Maciag, S. *et al.* (2022) Effects of freezing storage on the stability of maternal cellular and humoral immune components in porcine colostrum. [Vet Immunol Immunopathol. 254: 110520.](#)
39. Forner, R. *et al.* (2021) Distribution difference of colostrum-derived B and T cells subsets in gilts and sows. [PLoS One. 16 \(5\): e0249366.](#)
40. Muir, A. *et al.* (2024) Single-cell analysis reveals lasting immunological consequences of influenza infection and respiratory immunization in the pig lung. [PLoS Pathog. 20 \(7\): e1011910.](#)
41. Sedaghat-Rostami, E. *et al.* (2025) Porcine respiratory coronavirus as a model for acute respiratory disease: mechanisms of different infection outcomes. [J Immunol. 214 \(7\): 1643-1660.](#)
42. Jaudas, F. *et al.* (2025) Perinatal dysfunction of innate immunity in cystic fibrosis. [Sci Transl Med. 17 \(782\): eadk9145.](#)
43. Tng, P.Y.L. *et al.* (2025) Early disruption of the innate-adaptive immune axis in vivo after infection with virulent Georgia 2007/1 ASFV. [Discov Immunol. 4 \(1\): kyaf014.](#)

---

**Further Reading**

1. Piriou-Guzylack, L. (2008) Membrane markers of the immune cells in swine: an update. [Vet Res. 39: 54.](#)
2. Gerner W *et al.* (2015) Phenotypic and functional differentiation of porcine αβ T cells: current knowledge and available tools. [Mol Immunol. 66 \(1\): 3-13.](#)

---

**Storage**

This product is shipped at ambient temperature.  
Store at +4°C. DO NOT FREEZE.  
This product should be stored undiluted.

---

<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	This product is covered by U.S. Patent No. 10,150,841 and related U.S. and foreign counterparts
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20471 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA1971SBV570">https://www.bio-rad-antibodies.com/SDS/MCA1971SBV570</a>
<b>Regulatory</b>	For research purposes only

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M448544:260123'

**Printed on 20 Apr 2026**

---

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)