

Datasheet: MCA199B

BATCH NUMBER 163152

Description:	RAT ANTI MOUSE IgM HEAVY CHAIN:Biotin
Specificity:	IgM HEAVY CHAIN
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	LO-MM-9
Isotype:	IgG2a
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			250ng/ml
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on mouse IgM/Sepharose 4B from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃) 50% Glycerol
Approx. Protein Concentrations	IgG concentration 1 mg/ml

Immunogen	Polyclonal and monoclonal mouse IgM.
External Database Links	<p>UniProt:</p> <p>P01872 Related reagents</p> <p>P01873 Related reagents</p> <p>Entrez Gene:</p> <p>16019 Ighm Related reagents</p> <p>16019 Ighm Related reagents</p>
RRID	AB_321870
Fusion Partners	Spleen cells from immunized LOU/c rats were fused with cells of the IR983F rat myeloma cell line.
Specificity	Rat anti Mouse IgM, clone LO-MM-9 is directed against the mu chain of mouse IgM.
ELISA	This product may be used in an indirect ELISA or as a secondary reagent in a sandwich ELISA together with MCA1293 as the capture reagent and and Streptavidin:HRP (STAR5B) as the detection reagent.
References	<ol style="list-style-type: none"> 1. Bazin, H. <i>et al.</i> (1974) Three classes and four (sub)classes of rat immunoglobulins: IgM, IgA, IgE and IgG1, IgG2a, IgG2b, IgG2c. Eur J Immunol. 4 (1): 44-8. 2. Bazin, H. <i>et al.</i> (1984) Rat monoclonal antibodies. II. A rapid and efficient method of purification from ascitic fluid or serum. J Immunol Methods. 71: 9-16. 3. Busund, L.T. <i>et al.</i> (2003) Spontaneously formed tumorigenic hybrids of Meth A sarcoma cells and macrophages <i>in vivo</i>. Int J Cancer. 106 (2): 153-9. 4. Achtman AH <i>et al.</i> (2003) <i>Plasmodium chabaudi chabaudi</i> infection in mice induces strong B cell responses and striking but temporary changes in splenic cell distribution. J Immunol. 171 (1): 317-24. 5. Huang, C.H. <i>et al.</i> (2011) Airway inflammation and IgE production induced by dust mite allergen-specific memory/effector Th2 cell line can be effectively attenuated by IL-35. J Immunol. 187: 462-71. 6. Bermejo, D.A. <i>et al.</i> (2011) <i>Trypanosoma cruzi</i> infection induces a massive extrafollicular and follicular splenic B-cell response which is a high source of non-parasite-specific antibodies. Immunology. 132 (1): 123-33. 7. Sreekanthreddy, P. <i>et al.</i> (2015) A three-dimensional model of the human blood-brain barrier to analyse the transport of nanoparticles and astrocyte/endothelial interactions. F1000Res. 4: 1279.
Storage	<p>This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.</p>

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10328 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA199B>
10328

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M383379:210513'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)