

Datasheet: MCA2025T

**BATCH NUMBER 1702**

<b>Description:</b>	MOUSE ANTI HUMAN CD49b
<b>Specificity:</b>	CD49b
<b>Other names:</b>	INTEGRIN ALPHA 2 CHAIN, VLA-2
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	16B4
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	25 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			30ug/ml
Immunohistology - Frozen	▪			30ug/ml - 40ug/ml
Immunohistology - Paraffin (1)	▪			
ELISA	▪			10ug/ml
Immunoprecipitation	▪			25ug/ml - 50ug/ml
Western Blotting	▪			
Immunofluorescence	▪			
Functional Assays (2)	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

**(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.**

**Tris/EDTA pH9.0 is recommended for this purpose.**

**(2) This product contains sodium azide, removal by dialysis is recommended prior to use in functional assays. Dialysis cassettes [EQU003](#) are suitable for this purpose.**

<b>Target Species</b>	Human
-----------------------	-------

<b>Species Cross Reactivity</b>	Does not react with:Rat
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml
<b>Immunogen</b>	Purified human beta 1 preparation from HT1080 fibrosarcoma cell extract
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P17301</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">3673</a>    ITGA2    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	CD49B
<b>RRID</b>	AB_1102121
<b>Fusion Partners</b>	Spleen cells from immunized BALB/c mice were fused with cells of the X63/Ag8.653 mouse myeloma cell line
<b>Specificity</b>	<b>Mouse anti Human CD49b antibody, clone 16B4</b> recognizes human CD49b, also known as integrin $\alpha$ 2, collagen receptor, Platelet membrane glycoprotein Ia or VLA-2 subunit alpha. CD49b is a 1181 amino acid ~160 kDa single pass type-1 transmembrane glycoprotein possessing multiple <a href="#">FG-GAP</a> repeats and a single <a href="#">VWFA</a> domain. Mouse anti human CD49b antibody, clone 16B4 is identified as capable of immunoprecipitating a non-reducible alpha subunit from I <sup>125</sup> surface labeled human cell line extract. It is confirmed as specific to the alpha 2 subunit by relative expression of antigen on various cell lines by FACS, and its recognition of affinity purified alpha 2 beta 1 in dot blots.
<b>References</b>	<ol style="list-style-type: none"> <li>Fitter, S. <i>et al.</i> (1999) Transmembrane 4 superfamily protein CD151 (PETA-3) associates with beta 1 and alpha IIb beta 3 integrins in haemopoietic cell lines and modulates cell-cell adhesion. <a href="#">Biochem J. 338 ( Pt 1): 61-70.</a></li> <li>Sincock, P.M. <i>et al.</i> (1999) PETA-3/CD151, a member of the transmembrane 4 superfamily, is localised to the plasma membrane and endocytic system of endothelial cells, associates with multiple integrins and modulates cell function. <a href="#">J Cell Sci. 112 ( Pt 6):</a></li> </ol>

[833-44.](#)

3. Eaton, C.L. *et al.* (2010) Evaluation of the frequency of putative prostate cancer stem cells in primary and metastatic prostate cancer. [Prostate. 70 \(8\): 875-82.](#)
4. Gassmann, P. *et al.* (2009) CXCR4 regulates the early extravasation of metastatic tumor cells in vivo. [Neoplasia. 11: 651-61.](#)
5. Gassmann, P. *et al.* (2010) *In vivo* tumor cell adhesion in the pulmonary microvasculature is exclusively mediated by tumor cell--endothelial cell interaction. [BMC Cancer. 10: 177.](#)
6. Ivaska, J. *et al.* (2005) PKCepsilon-mediated phosphorylation of vimentin controls integrin recycling and motility. [EMBO J. 24: 3834-45.](#)
7. Pellinen, T. *et al.* (2006) Small GTPase Rab21 regulates cell adhesion and controls endosomal traffic of beta1-integrins. [J Cell Biol. 173: 767-80.](#)
8. Upla, P. *et al.* (2004) Clustering induces a lateral redistribution of alpha 2 beta 1 integrin from membrane rafts to caveolae and subsequent protein kinase C-dependent internalization. [Mol Biol Cell. 15: 625-36.](#)
9. Jokinen, J. *et al.* (2010) Molecular mechanism of alpha2beta1 integrin interaction with human echovirus 1. [EMBO J. 29: 196-208.](#)
10. Offtender, M. *et al.* (2003) Heregulin and retinoids synergistically induce branching morphogenesis of breast cancer cells cultivated in 3D collagen gels. [J Cell Physiol. 195: 260-75.](#)
11. Connors, W.L. and Heino, J. (2005) A duplexed microsphere-based cellular adhesion assay. [Anal Biochem. 337: 246-55.](#)
12. Connors, W.L. *et al.* (2007) Two synergistic activation mechanisms of alpha2beta1 integrin-mediated collagen binding. [J Biol Chem. 282: 14675-83.](#)
13. Vuoriluoto, K. *et al.* (2008) Syndecan-1 supports integrin alpha2beta1-mediated adhesion to collagen. [Exp Cell Res. 314: 3369-81.](#)
14. Holleran, B.J. *et al.* (2003) Differential recruitment of alpha2beta1 and alpha4beta1 integrins to lipid rafts in Jurkat T lymphocytes exposed to collagen type IV and fibronectin. [J Leukoc Biol. 73: 243-52.](#)
15. Rintanen, N. *et al.* (2012) Calpains promote  $\alpha 2 \beta 1$  integrin turnover in nonrecycling integrin pathway. [Mol Biol Cell. 23: 448-463.](#)
16. Pietiäinen, V. *et al.* (2004) Echovirus 1 endocytosis into caveosomes requires lipid rafts, dynamin II, and signaling events. [Mol Biol Cell. 15: 4911-25.](#)
17. Kankaanpää, P. *et al.* (2015) Cellular recognition and macropinocytosis-like internalization of nanoparticles targeted to integrin  $\alpha 2 \beta 1$ . [Nanoscale. 7 \(42\): 17889-901.](#)
18. Poulter, N.S. *et al.* (2017) Clustering of glycoprotein VI (GPVI) dimers upon adhesion to collagen as a mechanism to regulate GPVI signaling in platelets. [J Thromb Haemost. 15 \(3\): 549-564.](#)

---

**Storage**

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2025T">https://www.bio-rad-antibodies.com/SDS/MCA2025T</a> 10040
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight@488</a> , <a href="#">DyLight@550</a> , <a href="#">DyLight@650</a> , <a href="#">DyLight@680</a> , <a href="#">DyLight@800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M366066:200529'

Printed on 18 Jan 2024