

Datasheet: MCA2038F

BATCH NUMBER 1807

Description:	MOUSE ANTI CAT CD5:FITC
Specificity:	CD5
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	FE1.1B11
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Cat		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation	Purified IgG prepared by affinity chromatography on Protein A		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml		

RRID AB_323962

Specificity **Mouse anti Cat CD5 antibody, clone FE1.1B11** recognizes feline CD5, a 67kDa cell surface glycoprotein expressed by peripheral T lymphocytes and thymocytes.

Flow Cytometry Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Roccabianca, P. *et al.* (2000) Characterization of the diffuse mucosal associated lymphoid tissue of feline small intestine. [Vet Immunol Immunopathol. 75 \(1-2\): 27-42.](#)
2. Riondato, F. *et al.* (2003) Effects of interferon alpha (INF-alpha) therapy on peripheral blood lymphocyte subsets from FIV and FeLV naturally infected cats. [Vet Res Commun. 27: 429-32.](#)
3. Affolter, V.K. and Moore, P.F. (2006) Feline progressive histiocytosis. [Vet Pathol. 43: 646-55.](#)
4. Roccabianca, P. *et al.* (2006) Feline large granular lymphocyte (LGL) lymphoma with secondary leukemia: primary intestinal origin with predominance of a CD3/CD8(alpha) (alpha) phenotype. [Vet Pathol. 43: 15-28.](#)
5. Abkowitz, J.L. *et al.* (2009) In utero transplantation of monocytic cells in cats with alpha-mannosidosis. [Transplantation. 88: 323-9.](#)
6. Guzera, M. *et al.* (2016) The use of flow cytometry for immunophenotyping lymphoproliferative disorders in cats: a retrospective study of 19 cases. [Vet Comp Oncol. 14 Suppl 1: 40-51.](#)

Storage Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA2038F>
10041

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA928F\)](#)

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA1209F\)](#)

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

'M366082:200529'

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

Printed on 30 May 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)