

Datasheet: MCA2041A488

Description:	MOUSE ANTI BOVINE CD172a:Alexa Fluor® 488
Specificity:	CD172a
Other names:	SIRP ALPHA
Format:	ALEXA FLUOR® 488
Product Type:	Monoclonal Antibody
Clone:	CC149
Isotype:	IgG2b
Quantity:	100 TESTS/1ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/5 - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine		
Product Form	Purified IgG conjugated to Alexa Fluor®488 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Alexa Fluor®488	495	519
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.05 mg/ml		

Immunogen Bovine afferent veiled dendritic cells

External Database

Links

UniProt:

[O46631](#) [Related reagents](#)

Entrez Gene:

[327666](#) SIRPA [Related reagents](#)

Synonyms

MYD1, PTPNS1, SHPS1, SIRP

RRID

AB_2943340

Specificity

Mouse anti Bovine CD172a antibody, clone CC149 recognizes bovine CD172a, also known as MyD-1 antigen and SIRPA.

CD172a is a ~55 kDa single pass type 1 membrane protein belonging to the family of signal regulatory proteins (SIRP). CD172a has been identified as the receptor for CD47. Bovine CD172a is strongly expressed by splenic macrophages, monocytes and a subset of afferent lymph veiled cells (ALVC) and by dendritic cells in the skin.

Flow Cytometry

Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul

References

1. Smith, R. *et al.* (2003) A novel MyD-1 (SIRP-1alpha) signaling pathway that inhibits LPS-induced TNFalpha production by monocytes. [Blood.102: 2532-40.](#)
2. Brackenbury, L.S. *et al.* (2005) Identification of a cell population that produces alpha/beta interferon *in vitro* and *in vivo* in response to noncytopathic bovine viral diarrhea virus. [J Virol. 79: 7738-44.](#)
3. Price, S.J. & Hope, J.C.. (2009) Enhanced secretion of interferon-gamma by bovine gammadelta T cells induced by coculture with *Mycobacterium bovis*-infected dendritic cells: evidence for reciprocal activating signals. [Immunology. 126:201-8](#)
4. Waters, W.R. (2009) Signal regulatory protein alpha (SIRPalpha) cells in the adaptive response to ESAT-6/CFP-10 protein of tuberculous mycobacteria. [PLoS One. 4: e6414.](#)
5. Jensen, K. *et al.* (2014) Comparison of small interfering RNA (siRNA) delivery into bovine monocyte-derived macrophages by transfection and electroporation. [Vet Immunol Immunopathol. 158 : 224-32.](#)
6. Hussen J *et al.* (2014) The chemokine CCL5 induces selective migration of bovine classical monocytes and drives their differentiation into LPS-hyporesponsive macrophages *in vitro*. [Dev Comp Immunol. 47 \(2\): 169-77.](#)
7. Eger, M. *et al.* (2015) Impacts of parturition and body condition score on glucose uptake capacity of bovine monocyte subsets. [Vet Immunol Immunopathol. 166 \(1-2\): 33-42.](#)
8. Vachery N *et al.* (2015) An *in vitro* model to assess the immunosuppressive effect of tick saliva on the mobilization of inflammatory monocyte-derived cells. [Vet Res. 46 \(1\): 117.](#)
9. Tahoun, A. *et al.* (2015) Functional analysis of bovine TLR5 and association with IgA responses of cattle following systemic immunisation with H7 flagella. [Vet Res. 46: 9.](#)
10. Pridans, C. *et al.* (2016) A Csf1r-EGFP Transgene Provides a Novel Marker for Monocyte Subsets in Sheep. [J Immunol. 197 \(6\): 2297-305.](#)

11. Herry, V. *et al.* (2017) Local immunization impacts the response of dairy cows to *Escherichia coli* mastitis. [Sci Rep. 7 \(1\): 3441.](#)
12. Liu, J. *et al.* (2020) *Theileria annulata*. transformation altered cell surface molecules expression and endocytic function of monocyte-derived dendritic cells. [Ticks Tick Borne Dis. 11 \(3\): 101365.](#)
13. Okino, C.H. *et al.* (2020) A polymorphic CD4 epitope related to increased susceptibility to *Babesia bovis* in Canchim calves. [Vet Immunol Immunopathol. 230: 110132.](#)
14. Barut, G.T. *et al.* (2020) Transcriptomic profiling of bovine blood dendritic cells and monocytes following TLR stimulation. [Eur J Immunol. 50 \(11\): 1691-711.](#)
15. Kolar, Q.K. *et al.* (2020) Anatomical distribution of respiratory tract leukocyte cell subsets in neonatal calves. [Vet Immunol Immunopathol. 227: 110090.](#)
16. Tahoun, A. *et al.* (2020) Inflammasome Activation in Bovine Peripheral Blood-Derived Macrophages Is Associated with Actin Rearrangement. [Animals \(Basel\). 10 \(4\): 655.](#)
17. Ibeagha-Awemu, E.M. *et al.* (2021) Regionally Distinct Immune and Metabolic Transcriptional Responses in the Bovine Small Intestine and Draining Lymph Nodes During a Subclinical *Mycobacterium avium* subsp. *paratuberculosis* Infection. [Front Immunol. 12: 760931.](#)
18. Park, D.S. *et al.* (2021) Dynamic changes in blood immune cell composition and function in Holstein and Jersey steers in response to heat stress. [Cell Stress Chaperones. 26 \(4\): 705-20.](#)
19. Marzo, S. *et al.* (2022) Characterisation of dendritic cell frequency and phenotype in bovine afferent lymph reveals kinetic changes in costimulatory molecule expression. [Vet Immunol Immunopathol. 243: 110363.](#)
20. Casaro, S. *et al.* (2022) Flow cytometry panels for immunophenotyping dairy cattle peripheral blood leukocytes [Vet Immunol Immunopathol. 248: 110417.](#)
21. Cunha, P. *et al.* (2022) Simplified Approaches for the Production of Monocyte-Derived Dendritic Cells and Study of Antigen Presentation in Bovine. [Front Vet Sci. 9: 891893.](#)
22. Hong, S. *et al.* (2024) Impact of an Injectable Trace Mineral Supplement on the Immune Response and Outcome of *Mannheimia haemolytica* Infection in Feedlot Cattle. [Biol Trace Elem Res. Jun 10 \[Epub ahead of print\].](#)
23. Scatà, C.M. *et al.* (2024) Characterization of Cellular Immune System at Different Ages in Water Buffalo (*Bubalus bubalis*) [Journal of Buffalo Science. 13: 133-139.](#)
24. Muir, A. *et al.* (2024) Single-cell analysis reveals lasting immunological consequences of influenza infection and respiratory immunization in the pig lung. [PLoS Pathog. 20 \(7\): e1011910.](#)
25. Gilbert, F.B. *et al.* (2025) Expression of FcγR by bovine mononuclear blood leukocytes. [Dev Comp Immunol. 162: 105304.](#)
26. Cid de la Paz, M. *et al.* (2025) Exploring the impacts of colostrum on systemic immune development in dairy calves. [J Dairy Sci. 108 \(12\): 13890-13904.](#)
27. Grandoni, F. *et al.* (2023) Flow Cytometric Identification and Enumeration of Monocyte Subsets in Bovine and Water Buffalo Peripheral Blood. [Curr Protoc. 3 \(2\): e676.](#)

Further Reading

1. Howard, C.J. *et al.* (1999) Dendritic cells in cattle: phenotype and function. [Vet Immunol Immunopathol. 72 \(1-2\): 119-24.](#)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for

short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee 12 months from date of despatch

Acknowledgements This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchased product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA2041A488>

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL:Alexa Fluor® 488 \(MCA691A488\)](#)

Recommended Useful Reagents

[BOVINE DENDRITIC CELL GROWTH KIT \(PBP014KZZ\)](#)

[BOVINE DENDRITIC CELL GROWTH KIT \(PBP015KZZ\)](#)

[MOUSE ANTI BOVINE CD205:FITC \(MCA1651F\)](#)

[MOUSE ANTI BOVINE CD205 \(MCA1651GA\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M421968:230819'

Printed on 23 Apr 2026