

Datasheet: MCA2048

BATCH NUMBER 173991

Description:	MOUSE ANTI HUMAN CD222
Specificity:	CD222
Other names:	IGF-2 RECEPTOR
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	MEM-238
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	▪			1/25 - 1/50
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Membrane permeabilization is required for this application. The use of Leucoperm (Product Code [BUF09](#)) is recommended for this purpose.

Target Species	Human
Species Cross Reactivity	Reacts with: Rhesus Monkey N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from ascites

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% sodium azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant vaccinia virus containing CD222.
External Database Links	UniProt: P11717 Related reagents Entrez Gene: 3482 IGF2R Related reagents
Synonyms	MPRI
RRID	AB_323432
Specificity	Mouse anti Human CD222 antibody, clone MEM-238 recognizes human CD222, a 250 kDa transmembrane protein originally identified as the IGF II receptor. CD222 is ubiquitously expressed and is involved with internalization of a variety of ligands. Mouse anti Human CD222 antibody, clone MEM-238 recognizes an epitope located in the region between extracellular domains 2 and 5 (aa 192-697) of CD222 (Roberts <i>et al.</i> 2010). CD222 is primarily expressed intracellularly with a small percentage of molecules being located at the cell surface (5 - 10%).
Flow Cytometry	Use 10µl of the suggested working dilution to label 1x10 ⁶ cells in 100µl
References	- Godar, S. <i>et al.</i> (2002) CD222 (Mannose-6 phosphate / insulin-like growth factor II-receptor) Summary and Workshop Report. In Leucocyte Typing VII: White Cell Differentiation Antigens. Edited by Mason, D. <i>et al.</i> Oxford University Press. pp482-485. - Leksa, V. <i>et al.</i> (2002) The N terminus of mannose 6-phosphate/insulin-like growth factor 2 receptor in regulation of fibrinolysis and cell migration. J Biol Chem. 277 (43): 40575-82. - McCormick, P.J. <i>et al.</i> (2008) Palmitoylation controls recycling in lysosomal sorting and trafficking. Traffic. 9: 1984-97. - Rezgui, D. <i>et al.</i> (2009) Structure and function of the human Gly1619Arg polymorphism of M6P/IGF2R domain 11 implicated in IGF2 dependent growth. J Mol Endocrinol. 42: 341-56. - Roberts, R.C. <i>et al.</i> (2010) Mistargeting of SH3TC2 away from the recycling endosome causes Charcot-Marie-Tooth disease type 4C. Hum Mol Genet. 19: 1009-18. - Adachi, A. <i>et al.</i> (2010) Golgi-associated GSK3beta regulates the sorting process of

post-Golgi membrane trafficking. [J Cell Sci. 123: 3215-25.](#)

7. Osborne, D.G. *et al.* (2015) Monitoring receptor trafficking following retromer and WASH deregulation. [Methods Cell Biol. 130: 199-213.](#)

8. Hertel, A. *et al.* (2022) USP32-regulated LAMTOR1 ubiquitination impacts mTORC1 activation and autophagy induction. [Cell Rep. 41 \(10\): 111653.](#)

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
----------------	---

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA2048
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M412610:221114'

Printed on 21 Aug 2026