

Datasheet: MCA2050

Description:	MOUSE ANTI HUNTINGTIN
Specificity:	HUNTINGTIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	HDB4E10
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen (1)	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) Increased cytoplasmic staining, relative to nuclear, has been reported using formaldehyde as a fixative compared with acetone/methanol, see Wilkinson *et al.*

Species Cross Reactivity

Reacts with: Mouse, Human, Rabbit

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers	0.09% Sodium Azide
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant protein corresponding to amino acids 1844 - 2131 of huntingtin.
External Database Links	<p>UniProt: P42858 Related reagents</p> <p>Entrez Gene: 3064 HTT Related reagents</p>
Synonyms	HD, IT15
RRID	AB_322302
Specificity	<p>Mouse anti Huntingtin antibody, clone HDB4E10 reacts with an epitope corresponding to the HDB region (amino acids 1844 - 2131) of the huntingtin protein.</p> <p>Mouse anti Huntingtin antibody, clone HDB4E10 detects a ~350 kDa band on western blots but also detects smaller degradation products of huntingtin. The combined use of Mouse anti Huntingtin antibody, clones HDB4E10 (MCA2050) and HDC8A4 (MCA2051) demonstrate that huntingtin is enriched in neuronal cells in the brain (Jones 1999).</p>
Histology Positive Control Tissue	Brain
References	<ol style="list-style-type: none"> 1. Wilkinson, F. L. <i>et al.</i> (1999) Localization of rabbit huntingtin using a new panel of monoclonal antibodies. Molecular Brain Research. 69: 10-20. 2. Yao Y <i>et al.</i> (2015) A striatal-enriched intronic GPCR modulates huntingtin levels and toxicity. Elife. 4: 4. 3. Shirasaki, D.I. <i>et al.</i> (2012) Network organization of the huntingtin proteomic interactome in mammalian brain. Neuron. 75 (1): 41-57. 4. Jones, A.L. (1999) The localization and interactions of huntingtin. Philos Trans R Soc Lond B Biol Sci. 354 (1386): 1021-7. 5. Tao, T. & Tartakoff, A.M. (2001) Nuclear relocation of normal huntingtin. Traffic. 2 (6): 385-94. 6. Massai, L. <i>et al.</i> (2013) Development of an ELISA assay for the quantification of soluble huntingtin in human blood cells. BMC Biochem. 14: 34. 7. Ni, C.L. <i>et al.</i> (2016) Polyglutamine Tract Expansion Increases S-Nitrosylation of Huntingtin and Ataxin-1. PLoS One. 11 (9): e0163359.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for

short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR8...)	DyLight@800
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight@488 , DyLight@550 , DyLight@650 , DyLight@680 , DyLight@800 , FITC , HRP

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M383914:210513'

Printed on 17 Sep 2021