

Datasheet: MCA2164GA

BATCH NUMBER 1703

Description:	MOUSE ANTI CHICKEN CD4
Specificity:	CD4
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	2-35
Isotype:	IgG2b
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Chicken
Species Cross Reactivity	<p>Reacts with: Turkey</p> <p>N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Chicken embryonic thymocytes
Fusion Partners	Lymph node cells from immunised Balb/c mice were fused with cells of the SP2/0 myeloma cell line
Specificity	Mouse anti Chicken CD4, clone 2-35 recognizes the chicken homologue of human CD4, a ~64 kDa cell surface protein expressed by thymocytes and a subset of T cells (Luhtala et al. 1993). Mouse anti Chicken CD4, clone 2-35 has been demonstrated to recognize turkey CD4 (Li et al. 1998).
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul
References	<ol style="list-style-type: none"> Vainio, O. <i>et al.</i> (1989) Characterization of chicken CD4-expressing cells. Prog Clin Biol Res. 307: 45-56. Li, Z. <i>et al.</i> (1999) Cross-reactive anti-chicken CD4 and CD8 monoclonal antibodies suggest polymorphism of the turkey CD8alpha molecule. Poult Sci. 78 (11): 1526-31. Koskinen, R. <i>et al.</i> (1999) Cloning and modeling of the first nonmammalian CD4. J Immunol. 162 (7): 4115-21. Pavlova, S.P. <i>et al.</i> (2010) <i>In vitro</i> and <i>in vivo</i> characterization of glycoprotein C-deleted infectious laryngotracheitis virus. J Gen Virol. 91 (Pt 4): 847-57. Luhtala, M. <i>et al.</i> (1993) Analysis of chicken CD4 by monoclonal antibodies indicates evolutionary conservation between avian and mammalian species. Hybridoma. 12: 633-46. Rosa, A.C. <i>et al.</i> (2014) Isolation and molecular characterization of Brazilian turkey reovirus from immunosuppressed young poults. Arch Virol. 159 (6): 1453-7. Blohm, U. <i>et al.</i> (2016) Immunological Competence of Different Domestic Chicken Breeds Against Avian Influenza Infection. Avian Dis. 60 (1 Suppl): 262-8. Röhe, I. <i>et al.</i> (2017) Effect of feeding soybean meal and differently processed peas on the gut mucosal immune system of broilers. Poult Sci. 96 (7): 2064-73. Sachan, S. <i>et al.</i> (2015) Adjuvant potential of resiquimod with inactivated Newcastle disease vaccine and its mechanism of action in chicken. Vaccine. 33 (36): 4526-32. Kannan, T.A. <i>et al.</i> (2017) Age Related Changes in T Cell Subsets in Thymus and Spleen of Layer Chicken (<i>Gallus domesticus</i>) Int J Curr Microbiol App Sci. 6 (1): 15-19.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/MCA2164GA>
10040

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)

Goat Anti Mouse IgG (H/L) (STAR117...) [FITC](#)

Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M366272:200529'

Printed on 28 May 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)