

Datasheet: MCA2187F

Description:	RAT ANTI MOUSE CD22:FITC
Specificity:	CD22
Other names:	Lyb-8
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	OX-97
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

RRID AB_323529

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species Mouse

Product Form Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525

Preparation Purified IgG prepared by affinity chromatography on Protein G

Buffer Solution Phosphate buffered saline

Preservative 0.09% Sodium Azide
Stabilisers 1% Bovine Serum Albumin

Approx. Protein Concentrations IgG concentration 0.1 mg/ml

External Database Links
UniProt: [P35329](https://www.uniprot.org/entry/P35329) [Related reagents](#)

Entrez Gene:[12483](#) Cd22 [Related reagents](#)

Synonyms	Lyb-8, Siglec2
-----------------	----------------

Specificity	Rat anti Mouse CD22 antibody, clone OX-97 recognizes domain 2 of the murine cell surface glycoprotein CD22, which is a member of the sialoadhesin family. CD22 is expressed on a subpopulation of mature B-lymphocytes and may modulate signalling through the B-cell antigen receptor.
--------------------	--

Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul. The Fc region of monoclonal antibodies may bind non-specifically to cells expressing low affinity Fc receptors. This may be reduced by using SeroBlock FcR (BUF041A/B).
-----------------------	--

References	1. Duong, B.H. <i>et al.</i> (2010) Peripheral B cell tolerance and function in transgenic mice expressing an IgD superantigen. J Immunol. 184: 4143-58.
-------------------	--

Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	---

Guarantee	18 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[MOUSE SEROBLOCK FcR \(BUF041A\)](#)[MOUSE SEROBLOCK FcR \(BUF041B\)](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M344198:190111'

Printed on 20 May 2019