

Datasheet: MCA2226F

**BATCH NUMBER 167458**

<b>Description:</b>	MOUSE ANTI SHEEP MHC CLASS II DR MONOMORPHIC:FITC
<b>Specificity:</b>	MHC CLASS II DR MONOMORPHIC
<b>Format:</b>	FITC
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	37.68
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/5 - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Target Species

Sheep

### Species Cross Reactivity

Reacts with: Goat, Human, Bovine

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid

### Max Ex/Em

Fluorophore	Excitation Max (nm)	Emission Max (nm)
FITC	490	525

### Preparation

Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant

### Buffer Solution

Phosphate buffered saline

<b>Preservative</b>	0.09% sodium azide (NaN <sub>3</sub> )
<b>Stabilisers</b>	1% bovine serum albumin
<b>Approx. Protein Concentrations</b>	IgG concentration 0.1 mg/ml
<b>Immunogen</b>	Splenocytes from ATL mice.
<b>RRID</b>	AB_323966
<b>Specificity</b>	<p><b>Mouse anti Sheep MHC class II antibody, clone 37.68</b> recognizes a monomorphic epitope on ovine MHC class II DR molecules, constitutively expressed on antigen presenting cells such as dendritic cells, B lymphocytes, monocytes, macrophages, activated T lymphocytes and may be induced on a range of other cell types by interferon gamma.</p> <p>The major histocompatibility complex (MHC) is a cluster of genes some of which are important in the immune response to infections. In sheep, this complex is referred to as the ovine leukocyte antigen (OLA) region. There are 2 major types of MHC class IIa molecules encoded by the OLA which are DR and DQ each composed of an alpha and beta chain.</p>
<b>Flow Cytometry</b>	Use 10µl of the suggested working dilution to label 1 x 10 <sup>6</sup> cells in 100µl
<b>References</b>	<ol style="list-style-type: none"> <li>1. Puri, N.K. <i>et al.</i> (1987) Monoclonal antibodies to sheep MHC class I and class II molecules: biochemical characterization of three class I gene products and four distinct subpopulations of class II molecules. <a href="#">Vet Immunol Immunopathol. 15 (1-2): 59-86.</a></li> <li>2. Puri, N.K. &amp; Brandon, M.R. (1987) Sheep MHC class II molecules. II. Identification and characterization of four distinct subsets of sheep MHC class II molecules. <a href="#">Immunology. 62 (4): 575-80.</a></li> <li>3. Puri, N.K. <i>et al.</i> (1987) Monoclonal antibodies to sheep MHC class II molecules recognize all HLA-D or subsets of HLA-D region products. <a href="#">Hum Immunol. 20 (3): 195-207.</a></li> <li>4. Ballingall, K.T. <i>et al.</i> (1995) Analysis of the fine specificities of sheep major histocompatibility complex class II-specific monoclonal antibodies using mouse L-cell transfectants. <a href="#">Anim Genet. 26 (2): 79-84.</a></li> <li>5. Wang, Y. <i>et al.</i> (2017) Characterization of a secreted cystatin of the parasitic nematode <i>Haemonchus contortus</i> and its immune-modulatory effect on goat monocytes. <a href="#">Parasit Vectors. 10 (1): 425.</a></li> <li>6. Wang, Y. <i>et al.</i> (2020) Characterization of a rhodanese homologue from <i>Haemonchus contortus</i> and its immune-modulatory effects on goat immune cells <i>in vitro</i>. <a href="#">Parasit Vectors. 13 (1): 454.</a></li> <li>7. López-Fernández, A. <i>et al.</i> (2020) Effect of Allogeneic Cell-Based Tissue-Engineered Treatments in a Sheep Osteonecrosis Model. <a href="#">Tissue Eng Part A. 26 (17-18): 993-1004.</a></li> <li>8. Ehsan, M. <i>et al.</i> (2021) <i>Fasciola gigantica</i>. tegumental calcium-binding EF-hand protein 4 exerts immunomodulatory effects on goat monocytes. <a href="#">Parasit Vectors. 14 (1): 276.</a></li> </ol>
<b>Storage</b>	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for

short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10041 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2226F">https://www.bio-rad-antibodies.com/SDS/MCA2226F</a> 10041
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL:FITC \(MCA929F\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M413383:221122'

Printed on 15 Mar 2024