

Datasheet: MCA2243

Description:	MOUSE ANTI HUMAN KIR
Specificity:	KIR
Other names:	KILLER CELL IMMUNOGLOBIN-LIKE RECEPTORS
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	NKVFS1
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation	▪			
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide

Carrier Free Yes

Approx. Protein Concentrations IgG concentration 1.0 mg/ml

External Database Links

UniProt:

[P43626](#) [Related reagents](#)
[P43627](#) [Related reagents](#)
[P43628](#) [Related reagents](#)
[Q14954](#) [Related reagents](#)
[P43631](#) [Related reagents](#)
[P43632](#) [Related reagents](#)

Entrez Gene:

[3802](#) KIR2DL1 [Related reagents](#)
[3803](#) KIR2DL2 [Related reagents](#)
[3804](#) KIR2DL3 [Related reagents](#)
[3806](#) KIR2DS1 [Related reagents](#)
[3809](#) KIR2DS4 [Related reagents](#)
[100132285](#) KIR2DS2 [Related reagents](#)

Synonyms CD158A, CD158B1, CD158B2, CD158H, CD158I, CD158J, KIRCL23, KKA3, NKAT1, NKAT2, NKAT5, NKAT6, NKAT8

RRID AB_2130355

Specificity

Mouse anti Human KIR antibody, clone NKVFS1 recognizes KIR2D members of the killer cell immunoglobulin (Ig)-like receptor (KIR) family, CD158a, CD158b and P50.3. KIR2D family members are cell surface glycoproteins with two Ig domains, which are expressed on natural killer cells and some T cells.

Mouse anti Human KIR antibody, clone NKVFS1 recognizes the long and short forms CD158a and CD158b (KIR2DL, KIR2DS1 and KIR2DS2 respectively) and also p50.3 (KIR2DS4).

Mouse anti Human KIR antibody, clone NKVFS1 is reported to have functional activity, activating NK cell cytotoxicity via KIR2DS and inhibiting via KIR2DL forms..

Flow Cytometry Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Spaggiari, G.M. *et al.* (2002) Soluble HLA class I molecules induce natural killer cell apoptosis through the engagement of CD8: evidence for a negative regulation exerted by members of the inhibitory receptor superfamily. [Blood. 99 \(5\): 1706-14.](#)
2. Spaggiari, G.M. *et al.* (2002) Soluble HLA class I induces NK cell apoptosis upon the engagement of killer-activating HLA class I receptors through FasL-Fas interaction. [Blood. 100 \(12\): 4098-107.](#)

3. Older Aguilar, A.M. *et al.* (2010) Coevolution of killer cell Ig-like receptors with HLA-C to become the major variable regulators of human NK cells. [J Immunol. 185 \(7\): 4238-51.](#)
4. Patterson, S. *et al.* (2008) Human invariant NKT cells display alloreactivity instructed by invariant TCR-CD1d interaction and killer Ig receptors. [J Immunol. 181 \(5\): 3268-76.](#)
5. Hilton, H.G. *et al.* (2015) The production of KIR-Fc fusion proteins and their use in a multiplex HLA class I binding assay. [J Immunol Methods. 425: 79-87.](#)
6. Moesta, A.K. *et al.* (2009) Chimpanzees use more varied receptors and ligands than humans for inhibitory killer cell Ig-like receptor recognition of the MHC-C1 and MHC-C2 epitopes. [J Immunol. 182 \(6\): 3628-37.](#)
7. Van Der Ploeg, K. *et al.* (2017) Modulation of Human Leukocyte Antigen-C by Human Cytomegalovirus Stimulates KIR2DS1 Recognition by Natural Killer Cells. [Front Immunol. 8: 298.](#)
8. Wang, Y. *et al.* (2009) Characteristics of expanded CD4+CD28null T cells in patients with chronic hepatitis B. [Immunol Invest. 38 \(5\): 434-46.](#)
9. Wijaya, R.S. *et al.* (2020) Expansion of dysfunctional CD56-CD16+ NK cells in chronic hepatitis B patients. [Liver Int. Dec 23 \[Epub ahead of print\].](#)

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
	Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
--------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR8...)	DyLight@800
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight@488 , DyLight@550 , DyLight@650 , DyLight@680 , DyLight@800 , FITC , HRP

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL (MCA928)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M386028:210518'

Printed on 06 Jan 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)