

Datasheet: MCA2246

**BATCH NUMBER 153561**

<b>Description:</b>	MOUSE ANTI HUMAN AMH
<b>Specificity:</b>	AMH
<b>Other names:</b>	ANTI MULLERIAN HORMONE
<b>Format:</b>	Con S/N
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	5/6
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			1/20 - 1/40
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.**

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	<p>Reacts with: Mouse, Sheep, Squirrel monkey, Baboon</p> <p><b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.</p>

<b>Product Form</b>	Concentrated tissue culture supernatant - liquid
<b>Preparation</b>	Concentrated Tissue Culture Supernatant containing 0.1M Tris/HCl pH7.4 and 5-10% foetal calf serum.
<b>Preservative Stabilisers</b>	0.1% sodium azide (NaN <sub>3</sub> )
<b>Immunogen</b>	Synthetic peptide derived from human AMH (VPTAYAGKLLISLSEERISAHHPNMVATEC)
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P03971</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">268</a>    AMH    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	MIF
<b>RRID</b>	AB_2226471
<b>Fusion Partners</b>	Spleen cells from immunised T/O outbred mice were fused with cells of the SP2/0 myeloma cell line.
<b>Specificity</b>	<p><b>Mouse anti Human AMH, clone 5/6</b> recognizes human anti-mullerian hormone (AMH), originally classified as a foetal testicular hormone that inhibits Mullerian duct development. AMH is expressed post-natally by immature Sertoli cells, and to a lesser degree by granulosa cells. AMH plays a role in testicular differentiation and in the regulation of ovarian follicle growth.</p> <p>AMH is a member of the TGF beta superfamily. It is secreted as a homodimeric ~140 kDa disulphide linked precursor that is cleaved to release the mature ~30 kDa homodimer.</p>
<b>Histology Positive Control Tissue</b>	human ovary
<b>References</b>	<ol style="list-style-type: none"> <li>1. Van Saen, D. <i>et al.</i> (2010) Meiotic activity in orthotopic xenografts derived from human postpubertal testicular tissue. <a href="#">Hum Reprod. 26: 282-93.</a></li> <li>2. Gruijters, M.J <i>et al.</i> (2003) Anti-Müllerian hormone and its role in ovarian function. <a href="#">Mol Cell Endocrinol. 211 (1-2): 85-90.</a></li> <li>3. Weenen, C. <i>et al.</i> (2004) Anti-Mullerian hormone expression pattern in the human ovary: potential implications for initial and cyclic follicle recruitment. <a href="#">Mol Hum Reprod10: 77-83.</a></li> <li>4. Papanastasopoulos, P. <i>et al.</i> (2009) A case of complete androgen insensitivity syndrome presenting with incarcerated inguinal hernia: an immunohistochemical study. <a href="#">Fertil Steril. 92: 1169.e11-4.</a></li> <li>5. Campbell, B.K. (2009) The endocrine and local control of ovarian follicle development in the ewe <a href="#">Anim. Reprod. 6:159-71</a></li> <li>6. Walker, M.L. <i>et al.</i> (2009) Ovarian aging in squirrel monkeys (<i>Saimiri sciureus</i>).</li> </ol>

[Reproduction. 138: 793-9.](#)

7. Sobinoff, A.P. *et al.* (2011) Understanding the villain: DMBA induced pre-antral ovotoxicity involves selective follicular destruction and primordial follicle activation through PI3K/Akt and mTOR signalling. [Toxicol Sci. 123: 563-75.](#)
8. Van Saen, D. *et al.* (2011) Can pubertal boys with Klinefelter syndrome benefit from spermatogonial stem cell banking? [Hum Reprod. 27: 323-30.](#)
9. David, A. *et al.* (2012) Effect of cryopreservation and transplantation on the expression of kit ligand and anti-Müllerian hormone in human ovarian tissue [Hum Reprod. 27: 1088-95.](#)
10. Kevenaar, M.E. *et al.* (2006) Serum anti-mullerian hormone levels reflect the size of the primordial follicle pool in mice. [Endocrinology. 147: 3228-34.](#)
11. Campbell, B.K. *et al.* (2012) The role of anti-Müllerian hormone (AMH) during follicle development in a monovulatory species (sheep). [Endocrinology. 153: 4533-43.](#)
12. Amorim, C.A. *et al.* (2013) Successful vitrification and autografting of baboon (*Papio anubis*) ovarian tissue. [Hum Reprod. 28: 2146-56.](#)
13. Parlakgumus, H.A. *et al.* (2015) GnRH agonists and antagonists in rescue for cyclophosphamide-induced ovarian damage: friend or foe? [Arch Gynecol Obstet. 291 \(6\): 1403-10.](#)
14. ThemmenA, P.N. *et al.* (2016) The use of anti-Müllerian hormone as diagnostic for gonadectomy status in dogs. [Theriogenology. 86 \(6\): 1467-74.](#)
15. Saatcioglu, H.D. *et al.* (2016) Control of Oocyte Reawakening by Kit. [PLoS Genet. 12 \(8\): e1006215.](#)
16. Calvopina, J.H. *et al.* (2015) The Aorta-Gonad-Mesonephros Organ Culture Recapitulates 5hmC Reorganization and Replication-Dependent and Independent Loss of DNA Methylation in the Germline. [Stem Cells Dev. 24 \(13\): 1536-45.](#)
17. Camlin, N.J. *et al.* (2016) Maternal Smoke Exposure Impairs the Long-Term Fertility of Female Offspring in a Murine Model. [Biol Reprod. 94 \(2\): 39.](#)
18. Ohta, K. *et al.* (2012) Male differentiation of germ cells induced by embryonic age-specific Sertoli cells in mice. [Biol Reprod. 86 \(4\): 112.](#)
19. Bazzano, M.V. *et al.* (2015) Obesity induced by cafeteria diet disrupts fertility in the rat by affecting multiple ovarian targets. [Reprod Biomed Online. 31 \(5\): 655-67.](#)
20. Díaz, P.U. *et al.* (2018) Altered Expression of Anti-Müllerian Hormone during the Early Stage of Bovine Persistent Ovarian Follicles. [J Comp Pathol. 158: 22-31.](#)
21. Cacciottola, L. *et al.* (2020) Long-Term Advantages of Ovarian Reserve Maintenance and Follicle Development Using Adipose Tissue-Derived Stem Cells in Ovarian Tissue Transplantation. [J Clin Med. 9 \(9\) Sep 15 \[Epub ahead of print\].](#)
22. Lv, X. *et al.* (2020) Effects of Single Transplantation and Multiple Transplantation of Human Umbilical Cord Mesenchymal Stem Cells on the Recovery of Ovarian Function in the Treatment of Premature Ovarian Failure in Mice [Res Sq; Dec 01 \[Preprint, Epub ahead of print\].](#)
23. Wang, F. *et al.* (2013) Human amniotic epithelial cells can differentiate into granulosa cells and restore folliculogenesis in a mouse model of chemotherapy-induced premature ovarian failure. [Stem Cell Res Ther. 4 \(5\): 124.](#)
24. de Michele, F. *et al.* (2018) *In vitro* formation of the blood-testis barrier during long-term organotypic culture of human prepubertal tissue: comparison with a large cohort of pre/peripubertal boys. [Mol Hum Reprod. 24 \(5\): 271-82.](#)

**Storage** This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

<b>Guarantee</b>	12 months from date of despatch
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10451 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2246">https://www.bio-rad-antibodies.com/SDS/MCA2246</a> 10451
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> , <a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M413874:221128'

Printed on 24 Jun 2024