

Datasheet: MCA2270

BATCH NUMBER 160629

Description:	MOUSE ANTI HUMAN CD321
Specificity:	CD321
Other names:	JAM-1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	M.Ab.F11
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/10 - 1/50
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting	▪			1/100 - 1/1000
Immunofluorescence	▪			
Functional Assays (1)	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product contains sodium azide, removal by dialysis is recommended prior to use in functional assays. Dialysis cassettes [EQU003](#) are suitable for this purpose.

Target Species	Human
Product Form	Purified IgG - liquid
Buffer Solution	Phosphate buffered saline

Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Human platelet membranes.
External Database Links	<p>UniProt: Q9Y624 Related reagents</p> <p>Entrez Gene: 50848 F11R Related reagents</p>
Synonyms	JAM1, JCAM
RRID	AB_324422
Fusion Partners	Spleen cells from immunized Balb/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Specificity	<p>Mouse anti Human CD321 antibody, clone F11 recognizes human Junctional adhesion molecule A, also known as CD321, JAM-1, JAM-A, Platelet F11 receptor, Platelet adhesion molecule 1 or PAM-1. CD321 is a 272 amino acid ~35 kDa single pass type 1 transmembrane glycoprotein expressed by endothelium, epithelial tissues and some leukocytes (Ostermann et al. 2002). CD321 was originally identified as the F11 receptor expressed on platelets (Kornecki et al. 1990).</p> <p>Mouse anti Human CD321 antibody, clone F11 is able to activate platelets, leading to aggregation and granule secretion (Kornecki et al. 1990, Naik et al. 1995).</p>
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ml
Western Blotting	MCA2270 detects a band of approximately 35-40kDa in platelet cell lysates.
References	<ol style="list-style-type: none"> Sobocka, M.B. <i>et al.</i> (2000) Cloning of the human platelet F11 receptor: a cell adhesion molecule member of the immunoglobulin superfamily involved in platelet aggregation. Blood. 95 (8): 2600-9. Fraemohs, L. <i>et al.</i> (2004) The functional interaction of the beta 2 integrin lymphocyte function-associated antigen-1 with junctional adhesion molecule-A is mediated by the I domain. J Immunol. 173 (10): 6259-64. Immenschuh, S. <i>et al.</i> (2009) Transcriptional induction of junctional adhesion molecule-C gene expression in activated T cells. J Leukoc Biol. 85 (5): 796-803. Haarmann, A. <i>et al.</i> (2010) Evaluation of soluble junctional adhesion molecule-A as a biomarker of human brain endothelial barrier breakdown. PLoS One. 5 (10): e13568.
Storage	Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA2270 10040
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M366494:200529'

Printed on 02 Oct 2023