

Datasheet: MCA2271T

Description:	MOUSE ANTI HUMAN DESMOGLEIN 1
Specificity:	DESMOGLEIN 1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	27B2
Isotype:	IgG1
Quantity:	25 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA			▪	
Immunoprecipitation	▪			
Western Blotting	▪			1/100 - 1/1000
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant human desmoglein 1.
External Database Links	UniProt: Q02413 Related reagents
	Entrez Gene:

Synonyms	CDHF4
-----------------	-------

RRID	AB_2093298
-------------	------------

Fusion Partners	Spleen cells from immunised Balb/c mice were fused with cells of the mouse NS1 myeloma cell line.
------------------------	---

Specificity	<p>Mouse anti Human desmoglein 1 antibody, clone 27B2 recognizes human desmoglein 1, also known as Cadherin family member 4, Desmosomal glycoprotein 1, DG1 or Pemphigus foliaceus antigen. Desmoglein 1 is a 999 amino acid mature peptide with an an additional 26aa propeptide region and a 23 aa signal peptide. Desmoglein 1 is a ~165 kDa single pass type I transmembrane glycoprotein and core component of the desmosome cell-cell junction between epithelial cells.</p> <p>Desmoglein 1 is expressed in the suprabasal layer of the skin and is vital to epidermal integrity (Kottke <i>et al.</i> 2006, Biedermann <i>et al.</i> 2010).</p> <p>Mouse anti Human desmoglein 1 antibody, clone 27B2 does not recognize desmoglein-2 or desmoglein-3.</p>
--------------------	---

Histology Positive Control Tissue	Skin
--	------

Western Blotting	MCA2271 detects a band of approximately 165kDa in epithelial cell lysates.
-------------------------	--

References	<ol style="list-style-type: none">1. Wahl, J.K. 3rd (2002) Generation of monoclonal antibodies specific for desmoglein family members. Hybrid Hybridomics. 21 (1): 37-44.2. Mahoney, M.G. <i>et al.</i> (2010) Compound heterozygous desmoplakin mutations result in a phenotype with a combination of myocardial, skin, hair, and enamel abnormalities. J Invest Dermatol. 130 (4): 968-78.3. Wu, H. <i>et al.</i> (2000) Protection against pemphigus foliaceus by desmoglein 3 in neonates. N Engl J Med. 343 (1): 31-5.4. Otkarina, D.A. <i>et al.</i> (2011) IgG-induced clustering of desmogleins 1 and 3 in skin of patients with pemphigus fits with the desmoglein nonassembly depletion hypothesis. Br J Dermatol. 165 (3): 552-62.
-------------------	--

Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Rabbit Anti Mouse IgG (STAR8...) [DyLight@800](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@680](#),
[DyLight@800](#), [FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M366497:200529'

Printed on 11 Aug 2020

© 2020 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)