

Datasheet: MCA2336

BATCH NUMBER 173328

Description:	MOUSE ANTI DAZL
Specificity:	DAZL
Format:	S/N
Product Type:	Monoclonal Antibody
Clone:	3/11A
Isotype:	IgG1
Quantity:	2 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat, Cynomolgus monkey

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Tissue Culture Supernatant - liquid

Preparation	Tissue Culture Supernatant containing 0.1M Tris/HCl
Preservative Stabilisers	<0.1% Sodium Azide (NaN ₃)
Immunogen	Synthetic peptide corresponding to sequence within the C terminal domain of human DAZL (CRVHHFRRSRAMLKSV).
External Database Links	<p>UniProt: Q92904 Related reagents</p> <p>Entrez Gene: 1618 DAZL Related reagents</p>
Synonyms	DAZH, DAZL1, DAZLA, SPGYLA
RRID	AB_2292585
Fusion Partners	Spleen cells from immunised T/O outbred mice were fused with cells of the SP2/0 myeloma cell line.
Specificity	<p>Mouse anti Human DAZL antibody, clone 3/11A recognizes human Deleted in azoospermia-like, also known as DAZL, DAZ homolog, DAZ-like autosomal, Deleted in azoospermia-like 1 or SPGY-like-autosomal. DAZL is a 295 amino acid ~33 kDa member of the DAZ (deleted in azoospermia) family of RNA binding proteins. DAZL is expressed in fetal and adult testes and ovaries, and is believed to play a role in germ cell development. In adult germ cells, the expression of DAZL is predominantly localized to the cytoplasm.</p> <p>Mutations in this gene have been linked to severe spermatogenic failure and infertility in males (Lin et al. 2001)</p>
Histology Positive Control Tissue	Ovary or testis.
References	<ol style="list-style-type: none"> Forand, A. & Bernardino-Sgherri, J. (2009) A critical role of PUMA in maintenance of genomic integrity of murine spermatogonial stem cell precursors after genotoxic stress. Cell Res. 19: 1018-30. Elkin, N.D. (2010) Toxicant-induced leakage of germ cell-specific proteins from seminiferous tubules in the rat: relationship to blood-testis barrier integrity and prospects for biomonitoring. Toxicol Sci.117: 439-48. Barrios, F. et al. (2010) Opposing effects of retinoic acid and FGF9 on Nanos2 expression and meiotic entry of mouse germ cells. J Cell Sci. 123: 871-80. Forand, A. et al. (2009) Similarities and differences in the <i>in vivo</i> response of mouse neonatal gonocytes and spermatogonia to genotoxic stress. Biol Reprod. 80: 860-73. Anderson, R.A. et al. (2007) Conserved and divergent patterns of expression of DAZL, VASA and OCT4 in the germ cells of the human fetal ovary and testis. BMC Dev Biol. 7: 136. Aoki, T. and Takada, T. (2012) Bisphenol A modulates germ cell differentiation and

- retinoic acid signaling in mouse ES cells. [Reprod Toxicol. 34: 463-70.](#)
7. Yamauchi K *et al.* (2009) *In vitro* germ cell differentiation from cynomolgus monkey embryonic stem cells. [PLoS One. 4 \(4\): e5338.](#)
 8. Zogbi, C. *et al.* (2012) Gonocyte development in rats: proliferation, distribution and death revisited. [Histochem Cell Biol. 138 \(2\): 305-22.](#)
 9. Woods DC *et al.* (2013) Embryonic stem cell-derived granulosa cells participate in ovarian follicle formation *in vitro* and *in vivo*. [Reprod Sci. 20 \(5\): 524-35.](#)
 10. Rose CM *et al.* (2014) Dynamic changes in DNA modification states during late gestation male germ line development in the rat. [Epigenetics Chromatin. 7: 19.](#)
 11. Chen SR *et al.* (2013) Disruption of genital ridge development causes aberrant primordial germ cell proliferation but does not affect their directional migration. [BMC Biol. 11: 22.](#)
 12. Jobling MS *et al.* (2011) Effects of di(n-butyl) phthalate exposure on foetal rat germ-cell number and differentiation: identification of age-specific windows of vulnerability. [Int J Androl. 34 \(5 Pt 2\): e386-96.](#)
 13. Conrad, S. *et al.* (2014) Differential gene expression profiling of enriched human spermatogonia after short- and long-term culture. [Biomed Res Int. 2014: 138350.](#)
 14. Wang Y *et al.* (2015) Protein Arginine Methyltransferase 5 (Prmt5) Is Required for Germ Cell Survival During Mouse Embryonic Development. [Biol Reprod. pii: biolreprod.114.127308.](#)
 15. Bayne RA *et al.* (2015) GDF9 is Transiently Expressed in Oocytes before Follicle Formation in the Human Fetal Ovary and is Regulated by a Novel NOBOX Transcript. [PLoS One. 10 \(3\): e0119819.](#)
 16. Endo, T. *et al.* (2015) Periodic retinoic acid-STRA8 signaling intersects with periodic germ-cell competencies to regulate spermatogenesis. [Proc Natl Acad Sci U S A. 112 \(18\): E2347-56.](#)
 17. Tian-Zhong, M. *et al.* (2016) Critical role of Emx2 in the pluripotency - differentiation transition in male gonocytes via regulation of FGF9/NODAL pathway. [Reproduction. 151 \(6\): 673-81.](#)
 18. Chen, M. *et al.* (2019) Abnormal Meiosis Initiation in Germ Cell Caused by Aberrant Differentiation of Gonad Somatic Cell. [Oxid Med Cell Longev. 2019: 8030697.](#)
 19. Rocha-da-Silva, L. *et al.* (2019) Expression of genome defence protein members in proliferating and quiescent rat male germ cells and the Nuage dynamics. [PLoS One. 14 \(6\): e0217941.](#)
 20. Liang, J. *et al.* (2019) Induction of Sertoli-like cells from human fibroblasts by NR5A1 and GATA4. [Elife. 8:e48767.](#)
 21. Li, X. *et al.* (2025) Loss of full-length DAZL isoform disrupts PABPC1-dependent translational regulation and meiosis. [Cell Death Dis. 16 \(1\): 841.](#)

Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Health And Safety Information Material Safety Datasheet documentation #10053 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA2336>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Goat Anti Mouse IgG (STAR76...) [RPE](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@550](#),
[DyLight@650](#), [DyLight@680](#), [DyLight@800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M437786:250319'

Printed on 11 May 2026