

Datasheet: MCA2439GA

BATCH NUMBER 174488

Description:	MOUSE ANTI BOVINE IgG
Specificity:	IgG
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	IL-A2
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/2,500 - 1/15,000
Immunoprecipitation	▪			
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Purified bovine IgG
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the X63.Ag8.653 myeloma cell line.
Specificity	Mouse anti Bovine IgG, clone IL-A2 is specific for bovine Immunoglobulin G (IgG). The specificity of this antibody for bovine IgG1 and IgG2 has been demonstrated by ELISA.
References	1. Naessens, J. <i>et al.</i> (1988) Identification of isotypes and allotypes of bovine immunoglobulin M with monoclonal antibodies. Immunology. 63 (4): 569-74. 2. Buza, J. <i>et al.</i> (1997) CD5+ B lymphocytes are the main source of antibodies reactive with non-parasite antigens in <i>Trypanosoma congolense</i>- infected cattle. Immunology 92: 226-33. 3. VanDissel, J.T. <i>et al.</i> (2005) Bovine antibody-enriched whey to aid in the prevention of a relapse of <i>Clostridium difficile</i>-associated diarrhoea: preclinical and preliminary clinical data. J Med Microbiol. 54 (Pt 2): 197-205. 4. Dorneles, E.M. <i>et al.</i> (2015) Immune Response of Calves Vaccinated with <i>Brucella abortus</i> S19 or RB51 and Revaccinated with RB51. PLoS One. 10 (9): e0136696. 5. Tabynov, K. <i>et al.</i> (2016) Prime-booster vaccination of cattle with an influenza viral vector <i>Brucella abortus</i> vaccine induces a long-term protective immune response against <i>Brucella abortus</i> infection. Vaccine. 34 (4): 438-44. 6. Ryskeldinova, S. <i>et al.</i> (2021) Registered Influenza Viral Vector Based <i>Brucella abortus</i> Vaccine for Cattle in Kazakhstan: Age-Wise Safety and Efficacy Studies. Front Cell Infect Microbiol. 11: 669196. 7. Faria, A.R. <i>et al.</i> (2020) Immunoproteomics of <i>Brucella abortus</i> reveals potential of recombinant antigens for discriminating vaccinated from naturally infected cattle. Microb Pathog. 147: 104345. 8. Kelly, R.F. <i>et al.</i> (2022) Bovine Tuberculosis Antemortem Diagnostic Test Agreement and Disagreement in a Naturally Infected African Cattle Population. Front Vet Sci. 9: 877534.
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/MCA2439GA
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR12...)	RPE
Rabbit Anti Mouse IgG (STAR13...)	HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Recommended Useful Reagents

[PURIFIED BOVINE IgG \(PBP002\)](#)

[RABBIT F\(ab'\)₂ ANTI MOUSE IgG:HRP \(Human Adsorbed\) \(STAR13B\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M387156:210615'

Printed on 21 Aug 2026