

Datasheet: MCA2455A488

BATCH NUMBER 1605

Description:	MOUSE ANTI HUMAN C3a RECEPTOR:Alexa Fluor® 488
Specificity:	C3a RECEPTOR
Other names:	COMPLEMENT COMPONENT 3a RECEPTOR 1
Format:	ALEXA FLUOR® 488
Product Type:	Monoclonal Antibody
Clone:	hC3aRZ1
Isotype:	IgG1
Quantity:	100 TESTS/1ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Alexa Fluor®488 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Alexa Fluor®488	495	519
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin		
Approx. Protein	IgG concentration 0.05 mg/ml		

Concentrations

Immunogen Recombinant human C3aR protein.

External Database Links

UniProt:

[Q16581](#) [Related reagents](#)

Entrez Gene:

[719](#) C3AR1 [Related reagents](#)

Synonyms AZ3B, C3R1, HNFAG09

RRID AB_877451

Fusion Partners Cells from immunised BALB/c mice were fused with cells of the X63-Ag8.653 myeloma cell line.

Specificity **Mouse anti Human C3a Receptor antibody, clone hC3aRZ1** recognizes the receptor for human complement component C3a (C3aR), a ~55 kDa cell surface protein of the rhodopsin family of G protein-coupled seven-transmembrane receptors. C3a is an anaphylatoxin, which is a potent proinflammatory mediator that possesses immunoregulatory properties. C3a regulates inflammatory functions by binding to the C3a receptor. C3aR is expressed on various tissues and leucocytes including basophils, neutrophils, mast cells, eosinophils and activated T lymphocytes.

Expression levels of C3aR on human cells varies significantly between donors.

Flow Cytometry Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul.

References

1. Werfel, T. *et al.* (2000) Activated human T lymphocytes express a functional C3a receptor. [J Immunol. 165 \(11\): 6599-605.](#)
2. Soruri, A. *et al.* (2003) IL-4 down-regulates anaphylatoxin receptors in monocytes and dendritic cells and impairs anaphylatoxin-induced migration *in vivo*. [J Immunol. 170 \(6\): 3306-14.](#)
3. Zwirner, J. *et al.* (1999) Evaluation of C3a receptor expression on human leucocytes by the use of novel monoclonal antibodies. [Immunology. 97: 166-72.](#)
4. Kacani, L. *et al.* (2001) C5a and C5a(desArg) enhance the susceptibility of monocyte-derived macrophages to HIV infection. [J Immunol. 166: 3410-5.](#)
5. Raby, A.C. *et al.* (2009) Soluble TLR2 reduces inflammation without compromising bacterial clearance by disrupting TLR2 triggering. [J Immunol. 183: 506-17.](#)
6. Martel, C. *et al.* (2011) Requirements for membrane attack complex formation and anaphylatoxins binding to collagen-activated platelets. [PLoS One. 6: e18812.](#)
7. Hengartner, N.E. *et al.* (2015) Crucial Role of IL1beta and C3a in the *In Vitro*-Response of Multipotent Mesenchymal Stromal Cells to Inflammatory Mediators of Polytrauma. [PLoS One. 10: e0116772.](#)
8. Schreiber, A. *et al.* (2009) C5a receptor mediates neutrophil activation and ANCA-induced glomerulonephritis. [J Am Soc Nephrol. 20: 289-98.](#)

9. Török K (2012) Human T cell derived, cell-bound complement iC3b is integrally involved in T cell activation. [Immunol Lett. 143: 131-6.](#)
10. Patzelt, J. *et al.* (2015) Expression of anaphylatoxin receptors on platelets in patients with coronary heart disease. [Atherosclerosis. 238 \(2\): 289-95.](#)
11. Raby, A.C. *et al.* (2009) Soluble TLR2 reduces inflammation without compromising bacterial clearance by disrupting TLR2 triggering. [J Immunol. 183 \(1\): 506-17.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Acknowledgements This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA2455A488>
10041

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:Alexa Fluor® 488 \(MCA928A488\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M366929:200529'

Printed on 18 Jan 2024

