

Datasheet: MCA2468A647T

Description:	RAT ANTI HUMAN CMKLR1:Alexa Fluor® 647
Specificity:	CMKLR1
Other names:	CHEMOKINE-LIKE RECEPTOR 1
Format:	ALEXA FLUOR® 647
Product Type:	Monoclonal Antibody
Clone:	BZ332
Isotype:	IgG2a
Quantity:	25 TESTS/0.25ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Alexa Fluor® 647 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Alexa Fluor®647	650	665
Preparation	Purified IgG prepared by affinity chromatography on Protein G from supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	0.09% Sodium Azide		
	1%	Bovine Serum Albumin	
Approx. Protein Concentrations	IgG concentration 0.05 mg/ml		

Immunogen	Peptide corresponding to an epitope within the N-terminal region of human CMKLR1.
External Database Links	UniProt: Q99788 Related reagents Entrez Gene: 1240 CMKLR1 Related reagents
Synonyms	CHEMR23, DEZ
RRID	AB_2081374
Specificity	Rat anti Human CMKLR1 antibody, clone BZ332 recognizes human serpentine chemokine-like receptor 1 (CMKLR1), a novel chemoattractant receptor, which is also referred to Chem23 or DEZ. CMKLR1 is expressed on plasmacytoid dendritic cells (pDC) in blood but is not expressed on myeloid dendritic cells, lymphocytes, monocytes or granulocytes. Studies show that CMKLR1 is also expressed on dendritic cells, which have been generated <i>in vitro</i> from monocytes. The sole ligand for CMKLR1 is chemerin, a chemoattractant that exists at the mRNA level in a wide range of tissues, and as an inactive proform in serum and plasma. Studies suggest that chemerin, in plasma, is converted into a potent pDC attractant following blood coagulation and may be a process by which pDC's are recruited to sites of bleeding and tissue damage.
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	1. Zabel, B.A. <i>et al.</i> (2005) Chemokine-like receptor 1 expression and chemerin-directed chemotaxis distinguish plasmacytoid from myeloid dendritic cells in human blood. J Immunol. 174 (1): 244-51. 2. Maheshwari, A. <i>et al.</i> (2009) Epithelial cells in fetal intestine produce chemerin to recruit macrophages. Am J Physiol Gastrointest Liver Physiol. 297: G1-G10. 3. ChanY, L.T. <i>et al.</i> (2018) NK cells produce high levels of IL-10 early after allogeneic stem cell transplantation and suppress development of acute GVHD. Eur J Immunol. 48 (2): 316-29. 4. Li, L. <i>et al.</i> (2014) Chemerin-derived peptide C-20 suppressed gonadal steroidogenesis. Am J Reprod Immunol. 71 (3): 265-77.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	18 months from date of despatch.

Acknowledgements This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information Material Safety Datasheet documentation #10041 available at: 10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL:Alexa Fluor® 647 \(MCA6005A647\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M353907:190605'

Printed on 09 Feb 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)