

## Datasheet: MCA2474F

**BATCH NUMBER 1701**

<b>Description:</b>	RAT ANTI MOUSE CD71:FITC
<b>Specificity:</b>	CD71
<b>Other names:</b>	TRANSFERRIN RECEPTOR
<b>Format:</b>	FITC
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	8D3
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	0.1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse						
<b>Product Form</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid						
<b>Max Ex/Em</b>	<table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation Max (nm)</th> <th>Emission Max (nm)</th> </tr> </thead> <tbody> <tr> <td>FITC</td> <td>490</td> <td>525</td> </tr> </tbody> </table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	FITC	490	525
Fluorophore	Excitation Max (nm)	Emission Max (nm)					
FITC	490	525					
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant						
<b>Buffer Solution</b>	Phosphate buffered saline						
<b>Preservative</b>	0.09% Sodium Azide (NaN <sub>3</sub> )						
<b>Stabilisers</b>	1% Bovine Serum Albumin						
<b>Approx. Protein</b>	IgG concentration 0.1 mg/ml						

## Concentrations

---

**Immunogen** Mouse transformed endothelioma cell line t.end1.

---

## External Database Links

### UniProt:

[Q62351](#)    [Related reagents](#)

### Entrez Gene:

[22042](#)    Tfrc    [Related reagents](#)

---

**Synonyms** Tfrr

---

**Fusion Partners** Spleen cells from immunised Lewis rats were fused with the cells of the mouse NS0 myeloma cell line.

---

## Specificity

**Rat anti Mouse CD71 antibody, clone 8D3** recognizes mouse CD71, a 763 amino acid ~95 kDa single pass type II cell surface transmembrane glycoprotein, otherwise known as the transferrin receptor. CD71 is a major iron-binding protein, which plays a key role in the transport of iron into cells that require it. In mice, CD71 is widely expressed on a variety of cells, including Sertoli cells and cells which form the blood brain barrier (BBB) in the central nervous system.

Rat anti Mouse CD71 antibody, clone 8D3 recognizes native, soluble and denatured forms of murine CD71. Binding of the 8D3 antibody to CD71 does not inhibit the proliferation of cell lines tested, and does not interfere with the uptake of iron into cells.

Rat anti Mouse CD71 antibody, clone 8D3 has been used as a BBB transporter vector in mice and is suitable for studying CD71 expression, and iron uptake into different tissues, in the mouse ([Kissel et al. 1988](#)).

---

**Flow Cytometry** Use 10ul of the suggested working dilution to label  $10^6$

---

## References

1. Kissel, K. *et al.* (1998) Immunohistochemical localization of the murine transferrin receptor (TfR) on blood-tissue barriers using a novel anti-TfR monoclonal antibody. [Histochem cell Biol. 110: 63-72.](#)
2. Lee, H.J. *et al.* (2000) Targeting rat anti-mouse transferrin receptor monoclonal antibodies through blood-brain barrier in mouse. [J Pharmacol Exp Ther. 292 \(3\): 1048-52.](#)
3. Zhang, Y. *et al.* (2004) Intravenous RNA interference gene therapy targeting the human epidermal growth factor receptor prolongs survival in intracranial brain cancer. [Clin Cancer Res. 10 \(11\): 3667-77.](#)
4. Cabezón I *et al.* (2015) Trafficking of Gold Nanoparticles Coated with the 8D3 Anti-Transferrin Receptor Antibody at the Mouse Blood-Brain Barrier. [Mol Pharm. 12 \(11\): 4137-45.](#)
5. Lee, H.J. *et al.* (2002) Imaging gene expression in the brain *in vivo* in a transgenic mouse model of Huntington's disease with an antisense radiopharmaceutical and drug-targeting technology. [J Nucl Med. 43 \(7\): 948-56.](#)
6. Manich, G. *et al.* (2013) Study of the transcytosis of an anti-transferrin receptor

antibody with a Fab' cargo across the blood-brain barrier in mice. [Eur J Pharm Sci. 49 \(4\): 556-64.](#)

7. Sehlin, D. *et al.* (2016) Antibody-based PET imaging of amyloid beta in mouse models of Alzheimer's disease. [Nat Commun. 7: 10759.](#)

---

#### Further Reading

1. Jones AR & Shusta EV (2007) Blood-brain barrier transport of therapeutics via receptor-mediation. [Pharm Res. 24 \(9\): 1759-71.](#)
2. Pardridge WM (2007) Blood-brain barrier delivery of protein and non-viral gene therapeutics with molecular Trojan horses. [J Control Release. 122 \(3\): 345-8.](#)
3. Boado RJ *et al.* (2009) Engineering and expression of a chimeric transferrin receptor monoclonal antibody for blood-brain barrier delivery in the mouse. [Biotechnol Bioeng. 102 \(4\): 1251-8.](#)

---

#### Storage

Store at +4°C or at -20°C if preferred.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. This product is photosensitive and should be protected from light.  
Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

#### Guarantee

12 months from date of despatch

---

#### Health And Safety Information

Material Safety Datasheet documentation #10041 available at:  
<https://www.bio-rad-antibodies.com/SDS/MCA2474F>  
10041

---

#### Regulatory

For research purposes only

---

## Related Products

### Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL:FITC \(MCA1212F\)](#)

### Recommended Useful Reagents

[MOUSE SEROBLOCK FcR \(BUF041A\)](#)

[MOUSE SEROBLOCK FcR \(BUF041B\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M367008:200529'

**Printed on 25 Mar 2023**