

## Datasheet: MCA2477G

**BATCH NUMBER 163370**

<b>Description:</b>	MOUSE ANTI HUMAN IgG (Fc) (CH2 DOMAIN)
<b>Specificity:</b>	IgG (Fc) (CH2 DOMAIN)
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	8A4
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.2 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/50 - 1/100
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Cynomolgus monkey

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid

#### Preparation

Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide
<b>Carrier Free</b>	Yes
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0 mg/ml
<b>Immunogen</b>	Human IgG polyclonal.
<b>External Database Links</b>	<p><b>UniProt:</b></p> <p><a href="#">P01857</a>    <a href="#">Related reagents</a></p> <p><a href="#">P01859</a>    <a href="#">Related reagents</a></p> <p><a href="#">P01860</a>    <a href="#">Related reagents</a></p> <p><a href="#">P01861</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b></p> <p><a href="#">3500</a> IGHG1    <a href="#">Related reagents</a></p> <p><a href="#">3501</a> IGHG2    <a href="#">Related reagents</a></p> <p><a href="#">3502</a> IGHG3    <a href="#">Related reagents</a></p> <p><a href="#">3503</a> IGHG4    <a href="#">Related reagents</a></p>
<b>RRID</b>	AB_770116
<b>Fusion Partners</b>	Spleen cells from immunized Balb/c mouse were fused with cells of the P3-NS1/Ag1 myeloma cell line.
<b>Specificity</b>	<b>Mouse anti Human IgG (Fc) (CH2 Domain) antibody, clone 8A4</b> recognizes an epitope within the CH2 domain of the Fc region of all subclasses of human IgG.
<b>Flow Cytometry</b>	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells in 100ul.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Komatsuda, A. <i>et al.</i> (2008) Monoclonal immunoglobulin deposition disease associated with membranous features. <a href="#">Nephrol. Dial. Transplant. 23: 3888-94.</a></li> <li>2. Komatsuda, A. <i>et al.</i> (2013) Proliferative glomerulonephritis with discrete deposition of monoclonal immunoglobulin <math>\gamma</math>1 C(H) 2 heavy chain and <math>\kappa</math> light chain: A new variant of monoclonal immunoglobulin deposition disease. <a href="#">Pathol Int. 63: 63-7.</a></li> <li>3. Gehlsen, K. <i>et al.</i> (2012) Pharmacokinetics of engineered human monomeric and dimeric CH2 domains. <a href="#">MAbs. 4: 466-74.</a></li> <li>4. Kato, H. <i>et al.</i> (2015) Rapid Deterioration of the Renal Function Caused by the Coexistence of Intratubular Amyloidosis and Myeloma Cast Nephropathy. <a href="#">Intern Med. 54 (23): 3023-8.</a></li> </ol>
<b>Storage</b>	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for

short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2477G">https://www.bio-rad-antibodies.com/SDS/MCA2477G</a> 10040
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	<a href="#">RPE</a>
Goat Anti Mouse IgG IgA IgM (STAR87...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR76...)	<a href="#">RPE</a>
Rabbit Anti Mouse IgG (STAR13...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (STAR70...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (H/L) (STAR117...)	<a href="#">Alk. Phos.</a> , <a href="#">DyLight®488</a> , <a href="#">DyLight®550</a> , <a href="#">DyLight®650</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Rabbit Anti Mouse IgG (STAR9...)	<a href="#">FITC</a>
Goat Anti Mouse IgG (STAR77...)	<a href="#">HRP</a>
Goat Anti Mouse IgG (Fc) (STAR120...)	<a href="#">FITC</a> , <a href="#">HRP</a>

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M383804:210513'

Printed on 01 Feb 2024