

Datasheet: MCA2486

BATCH NUMBER 159595

Description:	MOUSE ANTI HUMAN MCP-1
Specificity:	MCP-1
Other names:	CCL2
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	2.2-4A4-1A11
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			10ug/ml
ELISA	▪			1.0 - 2.0ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.2 - 0.4ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 0.5 ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.

Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution
Immunogen	Recombinant human MCP-1/MCAF (PHP061)
External Database Links	<p>UniProt: P13500 Related reagents</p> <p>Entrez Gene: 6347 CCL2 Related reagents</p>
Synonyms	MCP1, SCYA2
RRID	AB_2071556
Fusion Partners	Spleen cells from immunized BALB/c X ICR F ₁ mice were fused with cells of a mouse myeloma cell line.
Specificity	<p>Mouse anti Human MCP-1 antibody, clone 2.2-4A4-1A11 recognizes human C-C motif chemokine 2, also known as Monocyte chemoattractant protein 1 (MCP-1), Monocyte chemotactic and activating factor, Monocyte secretory protein JE, Small-inducible cytokine A2 or CCL2, and does not cross react with MCP-2, MCP-3 or other structurally related chemokines. MCP-1 is a 99 amino acid, including a 23 aa signal peptide ~13 kDa chemotactic protein, a chemoattractant for monocytes and basophils but not neutrophils or eosinophils (Weber et al.1996).</p> <p>MCP-1 is produced by a number of cell types, including monocytes, macrophages, microglia, fibroblasts and epithelial cells (Cushing et al. 1990). MCP-1 influences T-cell immunity polarizing responses towards the Th2 phenotype (Karpus and Kennedy 1997). Variations in the CCL2 gene affects susceptibility to <i>Mycobacterium tuberculosis</i> infection (Flores-Villanueva et al. 2005).</p> <p>Mouse anti Human MCP-1 antibody, clone 2.2-4A4-1A11 recognizes both native and recombinant MCP-1 in conjunction with Rabbit anti Human MCP-1 polyclonal antibody in a sensitive sandwich ELISA assay.</p>
ELISA	This product may be used in an indirect ELISA or as a capture antibody in a sandwich ELISA together with AHP784B as the detection reagent and PHP061 as the standard.

Histology Positive Control Tissue	Human Breast Invasive Ductal Carcinoma
Western Blotting	This product may be used under either reducing or non-reducing conditions together with PHP061 as the positive control.
References	<ol style="list-style-type: none"> Johansson, U. <i>et al.</i> (2008) Formation of composite endothelial cell-mesenchymal stem cell islets: a novel approach to promote islet revascularization. Diabetes. 57: 2393-401. Nio-Kobayashi J <i>et al.</i> (2015) Regulated C-C motif ligand 2 (CCL2) in luteal cells contributes to macrophage infiltration into the human corpus luteum during luteolysis. Mol Hum Reprod. pii: gav028. Chuang, L.P. <i>et al.</i> (2016) Increased MCP-1 gene expression in monocytes of severe OSA patients and under intermittent hypoxia. Sleep Breath. 20 (1): 425-33.
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C.</p> <p>This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: https://www.bio-rad-antibodies.com/SDS/MCA2486 10162
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	HRP
Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®550 , DyLight®650 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M372572:200706'

Printed on 13 Feb 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)