

Datasheet: MCA2489

**BATCH NUMBER 159766**

<b>Description:</b>	MOUSE ANTI STREP-TAG CLASSIC
<b>Specificity:</b>	STREP-TAG CLASSIC
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	Strep-tag II
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			0.2 ug/ml - 0.5 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Synthetic Peptide
<b>Product Form</b>	Purified IgG - lyophilised
<b>Reconstitution</b>	Reconstitute with 0.2 ml water. The addition of 0.09% sodium azide is recommended for long term storage.
<b>Preparation</b>	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	None present

## Stabilisers

---

### Approx. Protein Concentrations

IgG concentration 0.5 mg/ml

---

### Immunogen

Strep-tag II peptide (WSHPQFEK) conjugated to Keyhole Limpet Hemocyanin.

---

### RRID

AB\_609795

---

### Specificity

**Mouse anti Strep-Tag Classic antibody, clone Strep-tag II**, also known as StrepMAB-Classic, recognizes Strep-tag II, a widely used tag in protein expression applications. This antibody recognizes both-C- and N-terminal Strep-tag II and is especially suited to Western blot applications.

---

### References

1. Renzi F *et al.* (2015) Glycan-Foraging Systems Reveal the Adaptation of *Capnocytophaga canimorsus* to the Dog Mouth. [MBio. 6 \(2\): .pii: e02507-14.](#)
2. Gordon, C.A. *et al.* (2015) NUSAP1 expression is upregulated by loss of RB1 in prostate cancer cells. [Prostate. 75 \(5\): 517-26.](#)
3. Mavrakis, M. *et al.* (2016) Purification of recombinant human and *Drosophila* septin hexamers for TIRF assays of actin-septin filament assembly. [Methods Cell Biol. 136: 199-220.](#)
4. Oda, S. *et al.* (2015) Crystal Structure of Marburg Virus VP40 Reveals a Broad, Basic Patch for Matrix Assembly and a Requirement of the N-Terminal Domain for Immunosuppression. [J Virol. 90 \(4\): 1839-48.](#)
5. Renzi, F. *et al.* (2015) Glycan-foraging systems reveal the adaptation of *Capnocytophaga canimorsus* to the dog mouth. [MBio. 6 \(2\): e02507.](#)

---

### Storage

Prior to reconstitution store at +4°C.  
After reconstitution store at -20°C.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted.  
Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

### Guarantee

12 months from date of despatch

---

### Acknowledgements

Sold under license of IBA

---

### Health And Safety Information

Material Safety Datasheet documentation #10294 available at:  
<https://www.bio-rad-antibodies.com/SDS/MCA2489>  
10294

---

### Regulatory

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...)

[RPE](#)

Goat Anti Mouse IgG IgA IgM (STAR87...) [HRP](#)  
Goat Anti Mouse IgG (STAR76...) [RPE](#)  
Goat Anti Mouse IgG (STAR70...) [FITC](#)  
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)  
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)  
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)  
Goat Anti Mouse IgG (STAR77...) [HRP](#)  
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®550](#),  
[DyLight®650](#), [DyLight®680](#), [DyLight®800](#),  
[FITC](#), [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M390673:210929'

**Printed on 07 Feb 2025**

---

© 2025 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)